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EU-Type Examination Certificate

For details on the EU-Type Examination Certificate, refer to the Hamilton Process Analytics website: www.hamiltoncompany.com
ATEX Certificate TÜV 03 ATEX 7005 X – Analog Sensors and Armatures

EU-Type Examination Certificate reference:

ATEX Certificate TÜV 03 ATEX 7005 X, Issue 11 – Analog Sensors and Armatures

Tab. 1 Hamilton sensor families

Polylite Pro
Polylite Pro VP
Polylite RX
Polylite Plus ***, XP ***
Polylite Plus *** VP, XP *** VP
Polylite Plus ORP
Polylite Plus RX XP VP
Polyplast Pro
Polyplast Pro RX
EasyFerm Plus ***
EasyFerm Plus *** VP, LEVP
EasyFerm Plus ORP
EasyFerm Bio ***
EasyFerm Bio *** VP
MecoTrode
MecoTrode VP
OxyTrode
FermoTrode
ChemoTrode
ChemoTrode VP
ChemoTrode P
ChemoTrode ORP
Liq-Glass PG

*** pH glass

Tab. 2a ATEX / IECEx approved armatures (not for ChemoTrodes or FermoTrodes)

	a-length 120	a-length 225
RetractoFit™ Series*	--	✓
RetractoFit™ Bio Series*	--	✓
FlexiFit™ Series*	✓	--
Retractex™ Series*	--	✓

* several part numbers possible

Tab. 2b ATEX / IECEx approved armatures for ChemoTrodes or FermoTrodes

	a-length 120	a-length 150	a-length 200 / 250
MasterFit™ 120 (Ref 237200)	✓	--	--
MasterFit™ 150 (Ref 237225)	--	✓	--
MasterFit™ 200 (Ref 237235)	--	--	✓

Tab. 3a Maximum process temperature for gas atmospheres

Type 1	P ₁ ≤ 50 mW	P ₁ ≤ 125 mW	P ₁ ≤ 230 mW	P ₁ ≤ 360 mW
T4	125 °C	119 °C	111 °C	100 °C
T5	90 °C	84 °C	76 °C	65 °C
T6	71 °C	58 °C	39 °C	15 °C

Type 2	P ₁ ≤ 50 mW	P ₁ ≤ 100 mW	P ₁ ≤ 150 mW	P ₁ ≤ 360 mW
T4	125 °C	121 °C	117 °C	100 °C
T5	90 °C	86 °C	82 °C	65 °C
T6	59 °C	41 °C	22 °C	Not allowed

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Operating Instructions for pH and ORP sensors for process applications

Every product for shipment or sent back for repair must be decontaminated. If working with hazardous liquids observe and carry out the maintenance procedures, paying particular attention to cleaning and decontamination. If the product becomes contaminated with biohazardous, radioactive or chemical material, it should be cleaned.

Safety instructions

This sensor is only to be used for the intended purpose and under safe conditions. Improper use or misuse can be dangerous. Assembly and maintenance should only be carried out by trained personnel. Since these sensors are made of glass, they should be handled with care. Sensors with a pressurized reference system (EasyFerm Plus, EasyFerm Bio) have a risk of spontaneous glass breakage. When working with this kind of sensor we recommend wearing protective glasses. When cleaning or regenerating the sensors with acids/ bases or solvents, protective glasses and gloves are recommended. Take care that the PG13.5 thread and the o-ring are not damaged while the sensor is screwed into the process. O-Rings are subject to wear and tear and should be replaced regularly (at least once per year).

If a high pressure sensor is operated for a considerable period of time above 10 bar, it should not be exposed to high temperatures (> 80° C) immediately afterwards. Air/gas bubbles may stick to the sensitive area of the sensor. As a consequence, the measurement value might be wrong. In case of a shutdown of the power supply (230 V) the measurement value could be wrong.

Preparing the sensor

Carefully remove the watering cap. Rinse the sensor with water. Check the interior of the pH glass membrane for air bubbles. Allow any bubbles to rise to the top by shaking the sensor gently. Refillable sensor: Prior to measuring, open the reservoir stopper; close it after measuring. Check level of electrolyte in the sensor.

If the sensor is going to be pressurized in an armature or housing, the pressure should be 0.5 bars above process pressure.

Mounting

pH and ORP sensors perform best when installed at a 15° or more angle from the horizontal, except the Polylite Plus VP sensors those can also be mounted upside down.

Screw the sensor into the housing with the PG thread. Ensure that the torque of 1.5 Nm will not be exceeded to avoid damaging of the o-ring.

Electrical wiring

The sensors are equipped with an SB, KB, or a VP connector head (see Tab. 4). Before connecting the sensor to the cable, check that the connections are clean and dry. Do not touch the electrical contacts! VP connectors especially should not be disconnected in moisture condensing environments. Unstable signals, low slope, or long response time could indicate a moist or contaminated connector. Clean the VP connector head with a paper towel moistened with ethanol. Dry the connector head after this procedure with a dry paper towel. If the sensor still has the described problems replace the VP cable.

Polylite Plus VP sensors are equipped with a solution ground. The solution ground is a platinum wire located behind the liquid junction within the reference electrode chamber. It is used to match the electronic ground of the transmitter with the potential of the solution and must be connected to the terminal «potential matching» (E+H, SIEMENS), «auxiliary electrode» (Knick, Mettler-Toledo), or «solution ground» (Yokogawa, Emerson). The solution ground should not be used to ground the solution. The current flowing through the solution ground should not exceed 1 µA. Many transmitters offer additional sensor diagnostics based on a solution ground contact. Polylite Plus VP sensors do not have a diaphragm. Therefore the signal of the diaphragm or reference resistance is irrelevant.

Calibration and measurement

For calibration always use unused Hamilton DuraCal pH or ORP buffers.

Storing the sensor

Sensors should be stored with the watering cap attached, containing 3M KCl solution (Ref 238036) or storage solution (Ref 238931). Sensors stored dry exhibit temporary drifting values. If the sensor dries out inadvertently, it can be placed in storage solution overnight to regenerate.

Cleaning

In general, you can use acids, alkaline solutions, and commonly used solvents for brief periods to clean electrodes with a glass shaft. Flush with water after cleaning. After cleaning, electrodes are likely to exhibit sluggish response times for a certain period, so place it in storage solution for 15 minutes after cleaning. Polyplast type electrodes: Note limited resistance to chemicals, see «General» below. Ceramic diaphragm (liquid junction): If protein contamination occurs, the electrode should be immersed in a solution of 0.4% HCl + 5 g/l pepsin for several hours. If blackening of the diaphragm is apparent (silver compounds), the electrode should be immersed in a solution of 0.4% HCl + 76 g/l thiourea. If working with hazardous liquids observe and carry out the maintenance procedures, paying particular attention to cleaning and decontamination.

Cleaning, assembly and maintenance should be performed by personnel trained in such work. Do not use any abrasive tissues or cleaning materials and do not use any cleaning chemicals other than described above. Before removing the sensor from the measuring setup, always make sure that the setup is pressure-less and cold and that no process medium can be accidentally spilled. When removing and cleaning the sensor, it is recommended to wear safety glasses and protective gloves.

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pH: immerse sensor for 10 min in 0.1 – 1M NaOH, then for 10 min in 0.1 – 1M HCl. After regeneration, place the sensor in storage solution for a further 15 min.

ORP: metal surfaces can be cleaned with slightly abrasive substances, such as toothpaste or very fine scouring powder.

Disposal

The design of Hamilton sensors minimizes environmental impact. According to the EU directive 2012/19/EU (and the local laws), sensors should be disposed of as electrical or electronic waste, and not as municipal waste. Alternatively, it can be sent back to Hamilton for disposal.

Sensors with pressurized reference systems (EasyFerm Plus, EasyFerm Bio), should be disposed in the supplied polystyrene packing. This will help prevent injury should the sensor break.

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General

The life cycle of sensors is determined by requirements regarding response time, zero point, and slope. Measurement conditions, especially high temperatures and aggressive measurement solutions, can shorten the life cycle. Under favorable conditions, the life cycle may be 1 – 3 years at room temperature and about 1 – 3 months at 90°C. There is also a slight aging factor during storage. Polyplast types have limited chemical resistance. Extended contact with acidic or alkaline solutions and alcohols should be avoided. Ethers, esters, ketones, or aromatic and halogenated hydrocarbons attack the shaft material and should therefore not be allowed to come in contact with the sensor.

Typ: High-quality cables increase the length of the life cycle and yield more stable measurement values. See also «Lab & Process Analytics» at www.hamiltoncompany.com

USE IN POTENTIALLY EXPLOSIVE ATMOSPHERES

ATEX / UKEX / IECEx marking

Gas: CE 0035 / UKEX 2571 Ⓢ II 1/2 G Ex Ia IIC T4/T5/T6 Ga/Gb
Staub (Polyplast only): CE 0035 / UKEX 2571 Ⓢ II 1/2 G Ex Ia IIB T4/T5/T6 Ga/Gb
Dust: CE 0035 / UKEX 2571 Ⓢ II 1/2 D Ex Ia IIC T x °C Da/Db

Manufacturer:	Hamilton Bonaduz AG, 7402 Bonaduz, Switzerland
EU-Type Examination Certificate:	TUV 03 ATEX 7005 X
UK Type Examination Certificate:	TUV 21 UKEX 7047 X
IECEx Certificate of Conformity:	IECEx TUR 14.0001 X

EU-Type Examination Certificate, UK Type Examination Certificate, and IECEx Certificate of Conformity can be downloaded from www.hamiltoncompany.com. Conditions described in the ATEX EU-Type Examination Certificate, UK Type Examination Certificate, or IECEx Certificate of Conformity must be adhered to. All ATEX specified pH and ORP electrodes see Section «EU-Type Examination Certificate».

ATTENTION! In case a gas atmosphere and a dust atmosphere are, or could be, present at the same time, the risk of explosion must be examined carefully, and special precautions may be necessary.

Assembly

- The operator of equipment in potentially explosive atmospheres is responsible for ensuring that all components of the system are certified for that area classification and are compatible with each other.
- Electrodes may only be used in armatures intended for their specific fitting length (a-length). For installation in gas atmospheres, the following armatures are free to use:
 - for all sensors from Tab. 1 with the exception of FermoTode and ChemoTode types: see Tab. 2a.
 - armatures for FermoTode and ChemoTode types: see Tab. 2b.As a general rule for installation in gas atmospheres, the glass surface in contact with the medium must be no more than 4 cm² for equipment group IIC, 25 cm² for equipment group IIB, and 50 cm² for equipment group IIA. For installations in dust atmospheres, special restrictions may apply when using armatures made of plastic material.
- ATTENTION:** Incorrect fitting or dismantling of the armature or the electrode may lead to unintended release of a potentially explosive atmosphere.
- O-Rings having a sealing function between Ex-Zone 0 and 1, 20 and 21 respectively must be replaced at each dismantling of the sensor.
- Sensors, wiring sections, transmitters and other required equipment are to be set up within a single potential equalizing system.
- When the sensor is assembled, the supplied ATEX / UKEX / IECEx sticker is to be attached to the cable in an easily visible place near the sensor. The sticker indicates that an ATEX / UKEX / IECEx approved sensor is mounted. The sticker must not be removed.

Connection to transmitter

The sensors listed in Tab. 1 are suitable for connection to an intrinsically safe electrical circuit with protection level Ia. The operator of the equipment must ensure that the allowable electrical values for the sensor all exceed those of the transmitter. None of the values cited for electrical power, voltage, and current may be exceeded in total (measuring and temperature circuit together). Electrical values for the sensors listed in table 1 are as follows: U = 24 V, I = 173 mA; The maximum process temperature for gas atmospheres is defined in Tab. 3a. For dust atmospheres, the maximum surface temperature of the sensor has to be observed. This surface temperature depends on the maximum electrical power of the transmitter, see Tab. 3b. The maximum allowable electric power P of the transmitter (measuring and temperature circuit together) depends on the sensor type (Tab. 1) and the desired process temperature.

Attention: The temperature limit given on the sensor must in any case be respected.

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