Hygienic Socket

HAMILTON 

Table of Contents

1	Introduction	4
2	Safety instructions	4
2.1	General safety instructions	4
2.2	Intended use	5
2.3	Basic safety	6
2.4	Other risks	6
3	Transportation and storage	7
4	Description of the Hygienic Socket	7
5	Assembly	8
5.1	General	8
5.2	Mounting the Hygienic Socket	8
5.3	Assembling the Armature	10
6	Operation	12
6.1	The «Life Guard» feature	12
7	Cleaning and maintenance	13
7.1	General cleaning information	13
7.2	Cleaning in place (CIP)	13
7.3	Cleaning the interior of the Hygienic Socket	13
7.4	Care of O-rings	15
7.5	Maintenance	15
8	Ordering information	16
8.1	Accessories	16
9	Technical Data	17

Hamilton Warranty

Please refer to the General Terms of Sales (GTS).

Important note

Copyright © 2021 Hamilton Bonaduz AG, Bonaduz Switzerland. All rights reserved.
The reproduction of any part of this document in any form is forbidden without the express written agreement of Hamilton Bonaduz AG.

Contents of this manual can be modified without previous announcement. Technical modifications reserved. Greatest possible care was used on the correctness of the information in this manual. If errors should be discovered nevertheless, Hamilton Bonaduz AG is pleased to be informed about it. Regardless of this, Hamilton Bonaduz AG cannot assume liability for any errors in this manual or for their consequences.

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1 Introduction

This manual is for the Hygienic Socket, developed especially for hygienic processes. The following versions of the Hygienic Socket are available:

Ref	Name
242 535	Hygienic Socket 1.4404
242 545	Hygienic Socket 1.4435
242 548	Hygienic Socket 1.4571
242 550	Hygienic Socket 2.4602

Hygienic Sockets are quality products, produced in accordance with the latest scientific and technical findings. Follow the instructions given here to ensure maximum safety and product longevity.

These instructions must be read, understood, and followed by all staff using the device. Hamilton assumes no responsibility for damage and/or operational disruption arising from failure to observe these instructions.

2 Safety instructions

2.1 General safety instructions

Hygienic Sockets and their parts are built using state-of-the-art technology, and are safe to operate. However, there can be risk to life and limb if users do not operate them correctly and appropriately.

A Hygienic Socket must only be used:

- For its intended purpose
- When in optimum condition

⚠ ATTENTION! Everyone involved in the operation and maintenance of a Hygienic Socket must carefully read and follow these operating instructions. Your personal safety may be at stake.

⚠ ATTENTION! The process pressure of 16 bar, and the process temperature of 140°C, must in no circumstances be exceeded.

⚠ ATTENTION! The Hygienic Socket must be maintained according to the instructions given in Section 7.5, «Maintenance».

The maximum process temperature permitted is limited not only by the Hygienic Socket itself, but also by the sensor installed in the Hygienic Socket. For this reason, always check the temperature specifications of the sensor to determine the maximum operating temperature.

Other relevant regulations must be complied with. These include, in particular:

- Regulations regarding safety and health protection in preparing materials and their use in the workplace
- Machine Directive
- Other national regulations
- Accident prevention regulations
- Any other relevant DIN and VDE directives

2.2 Intended use

Hygienic Sockets are designed for specific operating conditions (see Section 6) with respect to temperature, pressure, and mounting.

Intended use also includes conformity with the conditions set down by Hamilton for assembly, disassembly, initial operation, operation, and maintenance. Any other or further use is classified as unintended, and Hamilton accepts no responsibility for any ensuing damage.

2.3 Basic safety

Danger of explosion

Where maintenance, assembly, or service work is needed on the Hygienic Socket, ensure that no potentially explosive atmosphere arises.

Installation, initial operation, cleaning, maintenance, trouble shooting

- Prescribed maintenance and inspection tasks must be completed at the appropriate intervals, as given in Section 7.5, 'Maintenance'.
- Operators must be informed before the start of maintenance and servicing.
- Secure all system parts and media, connected to the Hygienic Socket against accidental process start-up.
- Sensors in the Hygienic Socket must be dismantled only in a de-pressurized condition.
- Defective parts must be replaced immediately. Use only original spare parts. There is no guarantee that parts from other manufacturers meet the operational and safety requirements of the Hygienic Socket.
- On completing maintenance, check that the inside of the Hygienic Socket is clean, and that none of the 'Life-Guard' boreholes are clogged or covered (such as by a cleaning adapter).

2.4 Other risks

⚠ ATTENTION! Even after all necessary safety measures have been complied with, there remains a risk of leaks or mechanical damage to the Hygienic Socket. In any place where there are seals or screws, there is the potential for gases or liquids to leak out undetectedly.


3 Transportation and storage

Always do the following:

- Check the Hygienic Socket package for completeness immediately after receipt.
- Notify in writing, without delay, any damage in transit or any missing parts.
- When returning a Hygienic Socket to a Hamilton partner, make sure it is in its original packaging.
- During storage protect the Hygienic Socket from moisture, direct heat sources (sun, heaters), contamination, and mechanical impact.

4 Description of the Hygienic Socket

The Hygienic Socket was developed especially for applications in biotechnology and pharmaceutical industry with stringent requirements for hygienic behavior.

5 Assembly

5.1 General


Electrical connections must be made only by qualified personnel.

Hamilton accepts no responsibility for damage resulting from incorrect assembly work by third parties.

If a Hygienic Socket is modified, Hamilton does not accept liability, and does not provide any further guarantee for the Hygienic Socket.

5.2 Mounting the Hygienic Socket

The steel sleeve must be welded into a circular cut-out of appropriate size (28 mm diameter) at its place of use. Installation depth can be 0 to 55 mm.

Hygienic Socket with sensor, mounted in wall.

The Hygienic Socket has optimized geometry between the sensor, O-ring, and internal sealing flange to provide a smooth, flush interface ideal for use in hygienic processes.

Too much heat from the welding process can deform the sealing flange and may lead to poor sealing of the sensor. Due to this potential issue, the following points need to be observed:

1. Only an experienced welding specialist should carry out the welding.
2. It is recommended to do a test weld to find out if the planned welding procedure is appropriate and the seal between Hygienic Socket and sensor is tight. If not, the welding procedure needs to be optimized and be approved again.
3. Flush welding of the hygienic Socket to the exposed metal surface is only possible with very thin plates (< 4 mm). For thicker metal plates it is recommended that the Hygienic Socket protrudes at least 25 mm past the surface.
4. Preparation of the welding groove (Double HV or K groove):
 - a. Hole diameter: 28.5 – 29.0 mm
 - b. Chamfer of plate: 30 – 45°
 - c. Thickness of root face around 1 mm


Figure 1

5. The heat input and arc energy of the welding process should be minimized to keep the temperature as low as possible to prevent deformation of the sealing flange. This can be achieved by a combination of the following procedures:
 - a. Cooling with forming gas (e.g. Argon)
 - b. Segmented welding (see figure 2):
 1. Welding to the root pass:
 - i. Welding of the first section (1)
 - ii. Welding of the opposite section (2)
 - iii. Cooling phase
 - iv. Welding of section 3
 - v. Welding of section 4

2. Welding of hot, fill and cover passes depending on plate thickness:
For every subsequent pass apply the same procedure as for the root pass
- c. Cooling with compressed air after each welding section.
- d. Welding with a welding current as small as possible, possibly by pulsed welding


Figure 2

NOTE: The CO₂NTROL sensor has been certified according to EHEDG Type EL Class I in combination with Hamilton Hygienic Socket REF 242535. Hence, it is important that in tanks/bioreactors, the cleaning device is positioned in such a way that the sensor is directly assessed and wetted for cleaning. The sensor should be installed flush with the process area. If mounting in a T-piece, the ratio between the depth of the upstand (L) and the diameter (D) of the upstand shall be $L/(D-d) < 1$ [d = 12mm diameter of sensor]. If welded adapters are used, the food contact surface must be smooth and the welding has to be done according to EHEDG Guideline 9 and 35 ($Ra < 0.8 \mu m$).

5.3 Assembling the Armature

NOTE: The Hygienic Socket is intended only for mounting sensors with a shaft length of 120 mm.


After the welded area has cooled down, do the following:

1. Ensure there is no damage to the sensor or the Hygienic Socket.

2. Check that all O-rings are in place in their appropriate grooves, and are free of damage.
3. Insert the seal pusher, and screw the sensor into the PG13,5 thread (locking torque see below).
4. Tightness of the installation must be verified.

When a sensor is screwed into the Hygienic Socket, the seal pusher presses the O-ring into the space between the steel sleeve and the sensor shaft. This seals the sensor in a manner that almost completely eliminates crevices and their associated risk of contamination. When the sensor is loosened, the O-ring relaxes again. This lowers the friction between the O-ring and the other parts of the assembly, and the process is no longer sealed. If this takes place while the process is still under pressure, the 'Life Guard' feature, explained below, with help to detect it.

Tip of a Hygienic Socket with loose sensor: the O-ring is not under pressure.

Tip of a Hygienic Socket with tightened sensor: the O-ring, under pressure, creates a seal with no crevices.

Recommended locking torque is 2 Nm (this is equivalent to a force of 10 N applied to the end of a 20 cm long spanner).

6 Operation

In EU countries, initial operation of the Hygienic Socket is permitted only when it has been verified that the Hygienic Socket conforms to the requirements of the EC directive on machines, the harmonized norms, the European norms, or the corresponding national norms.

Before initial operation of the Hygienic Socket:

- Ensure that the seal is tight.
- Ensure that all parts are in working order.

In addition to initial operation, perform the above checks:

- Immediately following any breakdown.
- When the device has been out of use for a long period.
- After repairs or maintenance.

⚠ ATTENTION! Contact with salt solutions and concentrated acids can result in corrosion.

Maximum process temperature: 140°C

Maximum process pressure: 16 bar

6.1 The «Life Guard» feature

The Hygienic Socket is equipped with a safety mechanism, the so-called 'Life Guard' boreholes, to prevent danger to personnel.

If the sensor is unscrewed while there is still pressure in the process, the O-ring between the sensor and the steel sleeve of the Hygienic Socket loses its seal. In this case, process medium enters the Hygienic Socket and exits through the 'Life Guard' boreholes (thereby alarming and protecting personnel) before the sensor is loosened entirely.

Fluid dripping from the cleaning/'Life-Guide' boreholes in the sleeve indicates that there is still pressure in the process. In this case, the sensor can be tightened again until such time as the process is no longer under pressure. Then the sensor can be safely removed.

Another advantage of the life guard feature is the control of the O-ring. If an O-ring is missing or broken during the process even though the sensor is mounted it will be shown because process medium will exit through the 'Life Guard' boreholes.

7 Cleaning and maintenance

Observe the general safety precautions given in these operating instructions when carrying out all maintenance, cleaning and servicing.

Both operational safety and the life expectancy of the Hygienic Socket are depending on appropriate maintenance and servicing, among other factors.

7.1 General cleaning information


The Hygienic Socket can be cleaned with standard cleaning agents. Never use scouring agents containing hard particles.

7.2 Cleaning in place (CIP)


If the installation into which the sensor is built is being cleaned, the parts of the Hygienic Socket and of the sensor that are normally in contact with the medium (now in contact with the cleaning liquid) can be cleaned in place, together with the rest of the installation.

7.3 Cleaning the interior of the Hygienic Socket

If the space between sensor and sleeve in the interior of the Hygienic Socket becomes contaminated, it can be cleaned easily by fitting cleaning adapters (not supplied by Hamilton) to the (non-threaded) boreholes in the wall of the steel sleeve, and injecting a suitable cleaning fluid. Alternatively, the Hygienic Socket can be steam sterilized.

Cleaning or steam-sterilizing the Hygienic Socket is possible by means of cleaning adapters inserted into the boreholes.


Dimensions of a borehole for cleaning adapters (dimensions shown in millimetres).

7.4 Care of O-rings

Check the O-rings regularly. When mounting a new O-ring, take care not to damage it, as a perfect seal can not be guaranteed with a damaged O-ring.

To replace the front O-ring, forcefully pull the sleeve from the Seal Pusher, then replace the O-ring and push the sleeve back onto the Seal Pusher.

Seal Pusher with O-ring sleeve removed. The O-ring can now be easily removed out of the sleeve.

After high temperature cycles (100 to 140°C) replace the O-ring each time the sensor is unscrewed. This is because the compression set of EPDM must be considered (DIN ISO 815, the deviation of the original O-ring form over time, due to high temperature and pressure).

Replace O-rings at least once a year.

7.5 Maintenance

Qualified personnel must draw up and implement a maintenance plan for the Hygienic Socket, according to the requirements of the installation.

Special attention must be given in the maintenance plan to the O-rings. O-rings used in Hygienic Sockets are consumable parts, and their life expectancy depends on process conditions. If the temperature frequently exceeds 100°C, the O-rings must be replaced each time the sensor is unscrewed, because they are easily deformed under conditions of high temperature. If replacement of O-rings according to the plan is not followed, leakage can occur.

- Check and replace O-rings as required; always replace at least once a year.
- Make sure the sensor is screwed in firmly the PG13,5 thread. (Weekly check recommended.)
- Perform CIP cleaning according to process requirements.

- Clean the Hygienic Socket by means of the ‘Life Guard’ boreholes, in case of contamination by the medium.
- Thoroughly clean the interior of the Hygienic Socket and the seal pusher after sensor breakage.
- Replace O-rings.


8 Ordering information

Hamilton bears no responsibility for damage arising from the use of non-original spare parts; this will void warranty.

8.1 Accessories

By default, Hygienic Sockets with Refs 242535, 242545 and 242548 are delivered with EPDM O-rings, Hygienic Socket 242550 is delivered with FFKM O-rings. O-rings of other materials are available on request.

DO Adapter (Ref 242538) is useful for mounting 120 mm oxygen or conductivity sensores, but only for sensors with a shaft diameter of 12 mm throughout the entire a-length. The DO Adapter may not be used for pH sensors or the Conducell UPW. Tightness of the installation after sensor mounting must be verified. The adapter is screwed onto the Hygienic Socket, and the sensor into the adapter. When in place, the sensor tip is almost flush with the Hygienic Socket, leading to excellent cleanability.


⚠ ATTENTION! The DO adapter must not be used with pH sensors or the Conducell UPW.

Seal Pusher Kits can be ordered separately (containing the Seal Pusher with EPDM O-rings which is put into the steel mantle).

Description	Ref
Hamilton DO Adapter	242538
Replacement Seal Pusher Kit for Hygienic Socket	242532
O-ring Set EPDM (10 O-rings of each size)	242595
O-ring Set FPM (10 O-rings of each size)	242596
O-ring Set SILICONE (10 O-rings of each size)	242597
O-ring Set FFKM (2 O-rings of each size)	242598

9 Technical Data

Ref 242535 / 242545 / 242548 / 242550	
Length	107 mm
Diameter	28 mm
This drawing is valid for Hygienic Socket of all Ref numbers. Hygienic Socket differ in the material of which the steel mantle is made, but not in their dimensions.	

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Copyright © 2021 Hamilton Bonaduz AG, Bonaduz Schweiz. Alle Rechte vorbehalten. Die Reproduktion irgendeines Teils dieses Dokuments in jeder beliebigen Form ist ohne die ausdrückliche schriftliche Zustimmung der Hamilton Bonaduz AG untersagt.
© 2021 Hamilton Bonaduz AG. All rights reserved.