OXYFERM FDA ARC

OXYGOLD G ARC

OXYGOLD B ARC

Operating Instructions

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Important note

Copyright © 2011 HAMILTON Bonaduz AG, Bonaduz Switzerland. All rights reserved. The reproduction of any part of this document in any form is forbidden without the express written agreement of HAMILTON Bonaduz AG.

Contents of this manual can be modified without previous announcement. Technical modifications reserved. Greatest possible care was used on the correctness of the information in this manual. If errors should be discovered nevertheless, HAMILTON Bonaduz AG is pleased to be informed about it. Regardless of this, HAMILTON Bonaduz AG cannot assume a liability for any errors in this manual or for their consequences.

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| | |
|---|----|
| IMPORTANT NOTE | 2 |
| WICHTIGER HINWEIS | 2 |
| 1. INTRODUCTION | 4 |
| 2. INTENDED USE..... | 4 |
| 3. SAFETY INSTRUCTIONS | 5 |
| 4. INITIAL OPERATION..... | 5 |
| 5. ELECTRICAL CONNECTION | 6 |
| ELECTRICAL CONNECTION OF THE 4-20 mA CURRENT INTERFACE..... | 7 |
| EXAMPLES OF CIRCUIT ARRANGEMENT | 8 |
| ELECTRICAL CONNECTION FOR THE DIGITAL RS485 INTERFACE | 9 |
| EXAMPLES OF CIRCUIT ARRANGEMENT | 10 |
| 6. CONFIGURATION AND MONITORING OF THE SENSOR..... | 11 |
| 7. PREPARATION FOR THE MEASUREMENT..... | 11 |
| 8. REMOVAL OF THE SENSOR..... | 12 |
| 9. STERILIZATION, AUTOCLAVING, CIP PURIFICATIONS..... | 12 |
| 10. TESTING AND MAINTENANCE | 13 |
| MAINTENANCE | 14 |
| SENSOR CHECK AND MAINTENANCE..... | 14 |
| REPLACING THE ELECTROLYTE AND THE MEMBRANE CARTRIDGE | 15 |
| REPLACING THE CATHODE | 16 |
| SELF-DIAGNOSIS FUNCTIONS..... | 17 |
| 11. DISPOSAL..... | 18 |
| 12. PARTS AND ACCESSORIES..... | 18 |
| 13. TECHNICAL DATA | 20 |

Operating Instructions OXYFERM FDA ARC and OXYGOLD ARC

1. Introduction

This manual refers to the OXYFERM FDA ARC and OXYGOLD ARC range of electrochemical, sterilizable oxygen sensors from HAMILTON Bonaduz AG. These sensors are compatible with all other components of the HAMILTON ARC System, including a complete family of intelligent sensors for pH, dissolved oxygen and conductivity measurements in a process control.

| Designation | Order number P/N |
|-----------------------|------------------|
| OXYFERM FDA ARC 125 | 243 100 |
| OXYFERM FDA ARC 160 | 243 101 |
| OXYFERM FDA ARC 225 * | 243 102 |
| OXYFERM FDA ARC 325 | 243 103 |
| OXYFERM FDA ARC 425 | 243 104 |
| OXYFERM FDA ARC XL | 243 140-OP** |
| OXYGOLD G ARC 120 | 243 110 |
| OXYGOLD G ARC 225* | 243 111 |
| OXYGOLD B ARC 120 | 243 115 |
| OXYGOLD B ARC 225* | 243 116 |

* The OXYFERM FDA ARC and OXYGOLD ARC 225 have an actual shaft length of 215 mm in order to ensure a good rinsing in replaceable armatures such as RETRACTEX.

** O-ring position.

HAMILTON sensors are quality products which are manufactured according to the latest knowledge. Follow the instructions given here and you can be sure of maximum safety and durability.

These instructions should be read, understood and followed by all staff using the device. HAMILTON can assume no responsibility for damage and operational disruptions arising from failure to observe these instructions.

2. Intended use

OXYFERM FDA ARC and OXYGOLD ARC sensors were developed for the measurement of the partial pressure of oxygen in liquids as well as the following derived measurement parameters:

- Percent volume of oxygen in air
- Percent air saturation by oxygen
- Concentration of oxygen in liquids

Main application for OXYFERM FDA ARC sensors is the biotechnology and pharmaceutical industry. For it, the sensor exhibits an unusual long-term stability even at frequent sterilization or autoclaving.

The OXYGOLD G ARC sensor is optimized for the measurement of trace quantities of dissolved oxygen. Its lower detection limit is 1 ppb. Typical application is the measurement of boiler feed water.

The OXYGOLD B ARC sensor is designed to measure the trace quantities of oxygen in conditions of high CO₂ partial pressure (beverage industry).

3. Safety instructions

OXYFERM FDA ARC and OXYGOLD ARC sensors must only be used for the intended application and in optimum safety and running conditions. The specifications such as temperature or pressure being defined in the chapter 'Technical data' may not be exceeded under any circumstances. Threats are imminent if the sensor is not operated correctly or appropriately.

Assembly and maintenance must be performed only by trained personnel.

Be careful that the process connection and O-ring are not damaged during the screwing in into the process. O-rings are wearing parts which regularly have to be exchanged, at the latest after one year. Even where all necessary safety measures have been complied with, there are still further risks involving leaks or mechanical damage to the armature. Where there are seals or screws, gases or liquids can leak out undetected.

Before removing the sensor from its measurement setup, always make sure that no process medium can be accidentally spilled.

The built-in temperature sensor can only be used for monitoring the sensor conditions, not for controlling the process temperature.

When changing the electrolyte, it is recommended that you wear safety goggles and protective gloves. Contact of the electrolyte with skin, eyes or mucous membrane is to be avoided (see Safety Data Sheet p/n 608914 for OXYFERM ARC and OXYGOLD G ARC type sensors and Safety Data Sheet p/n 608929 for OXYGOLD B type sensors).

4. Liability

The liability of HAMILTON Bonaduz AG is detailed in the document 'General Terms and Conditions of Sale and Delivery (GTS)', chapter 12.

Hamilton is expressly not liable for direct or indirect losses arising from use of the sensors. It must in particular be insured in this conjunction that malfunctions can occur on account of the inherently limited useful life of sensors contingent upon their relevant applications. The user is responsible for the calibration, maintenance and regular replacement of the sensors. In the case of critical sensor applications, Hamilton recommends using back-up measuring points in order to avoid consequential damages. The user is responsible for taking suitable precautions in the event of a sensor failure.

5. Initial operation

Control the OXYFERM FDA or OXYGOLD ARC sensor during unpacking according to possible defects. Rejected sensors are to be sent to your HAMILTON merchant in the original packaging.

Fig. 1: Components of electrochemical ARC sensor with 12 mm shaft:

1: Socket head with VP8 connector, 2: ARC logo, serial and part number, 3: PG 13.5 thread, 4: steel shaft with heat number, 5: membrane cap, 6: silicone membrane

OXYFERM FDA and ARC OXYGOLD ARC sensors are delivered directly from the factory already pre-calibrated (according to the operational specifications). The integrated 4–20 mA analog interface and RS485 digital interface (Modbus RTU) are configured according to factory defaults. You can find full details, including serial number and most important specifications, on the certificate provided with each sensor.

ATTENTION !

OXYFERM FDA ARC and OXYGOLD ARC sensors are not intended for use with any transmitter. They can send digital or 4-20 mA signals directly.

Avoid electrical damages at the sensor! Follow the instructions in the chapters ‘Electrical Connection’ and ‘Configuration of the sensor’ in order to adjust the sensor to your application and desires!

Before the assignment of the sensor for measurement, control or regulation you should examine the sensor configuration by means of an operational test.

6. Electrical connection

OXYFERM FDA ARC and OXYGOLD ARC sensors are fitted with a VP8 socket head. The eight golden contacts are denoted as pin A, pin B... pin H. For easy identification of each pin the head has a mark between pin A and pin B.

Fig. 2: VP pins

For connecting OXYFERM FDA ARC and OXYGOLD ARC sensors, always use HAMILTON VP 8.0 cables. These are available in a range of in different lengths.

| VP pin | Function |
|--------|-------------------------------------|
| A | 4-20 mA interface (mA interface #2) |
| B | 4-20 mA interface (mA interface #1) |
| C | Power supply: +24 VDC (7 to 30 VDC) |
| D | Power supply: ground |
| G | RS-485 (A) |
| H | RS-485 (B) |

Table 1: VP pin layout

Electrical connection of the 4-20 mA current interface

The 4–20 mA interface enables direct connection of the OXYFERM FDA ARC and OXYGOLD ARC sensor to a data recorder, indicator, control unit or PLC with analog I/O. Apart from the two wires used for 4–20 mA current loop no additional equipment is required for analog signaling.

When using the 4–20 mA interface, pins have the following designations with respect to VP cable conductor colors:

| OXYFERM FDA ARC, OXYGOLD ARC | VP pin | VP8 cable |
|---|--------|---------------------------------------|
| 4–20 mA three-wire interface, functions as a current sink. It regulates the input current according to the sensor measurements. This interface needs a separate power supply (pin B or C). Factory default is assigned to a temperature measurement. | A | Coaxial core black transparent. |
| 4–20 mA two-wire interface, functions as a current sink. It regulates the input current according to the sensor measurements. This interface can be powered directly from 4–20 mA two-wire current loop. Nominal power of 60 mW must be provided. Factory default is assigned to a oxygen mesurement. | B | Coaxial shield black. |
| Power supply: +24 VDC (7 to 30 VDC). Max. power consumption: 150 mW. | C | Coaxial core red transparent. |
| Power supply: Ground. | D | Coaxial shield red. |

Table 2: 4-20 mA interface pin layout

The 4–20 mA interface is configured at the factory with the range, value and unit for the measurements as indicated on the certificate. Follow the instructions in the section entitled 'Configuration and monitoring of the sensor' to adjust the sensor according to the requirements of your application.

Examples of circuit arrangement


Fig. 3: Two-wire loop wiring diagram for the 4–20 mA interface (mA interface #1). In this wiring scheme the power is not supplied to the sensor VP pin C, wherefore the wiring is not applicable to a sensor with ARC Wi Sensor Adapter.


Fig. 4: Three-wire loop wiring diagram for the 4–20 mA interfaces. The figure represents both 4–20 mA interfaces at pin A and pin B.


Fig. 5: The safest form of wiring, using an isolation amplifier. The figure represents both 4–20 mA interfaces at pin A and pin B. (For detailed technical advice, please contact the technical support at HAMILTON.)

Electrical connection for the digital RS485 interface

The digital RS485 interface enables communication with OXYFERM FDA ARC and OXYGOLD ARC sensor for performing measurements, for calibrating the sensor and for changing the sensor's configuration parameters. OXYFERM FDA ARC and OXYGOLD ARC sensors are always connected to digital controlling devices as a Modbus slave. To function, they require a power supply (VP8 pins C and D, see below). The section entitled 'Configuration and monitoring of the sensor' describes operation in digital mode.

By using the correct access password the ARC System operator can adapt ARC sensors to many tasks by:

- Selecting the 4–20 mA interface.
- Scaling (configuring) the 4–20 mA interface.
- Selecting the measured parameter:
 - Dissolved oxygen: partial pressure of oxygen (mbar); % oxygen air saturation (%-sat); volume % of oxygen (%-vol); concentration in liquid (µg/l, mg/l, ppb, ppm).
 - Temperature: °C, °F, K.

In addition, operators can read sensor information from the RS485 interface such as:

- The sensor's serial number, part number (P/N) and manufacturing number (WO).
- The sensor's firmware version.
- The sensor's status (e.g. operation hours, number of cleanings and sterilizations, warnings and errors).

Additional information:

The Modbus RTU communication protocol corresponds to the Modbus-IDA standard (see <http://www.modbus.org>). OXYFERM FDA ARC and OXYGOLD ARC sensors use an open register set developed by HAMILTON. Additional information about the register content and structure can be found at <http://www.hamiltoncompany.com>.

ATTENTION!

Because all sensors are delivered with factory-default settings, each sensor must be configured for its specific application before first use (see the section entitled 'Configuration of the sensor').

The pins for the digital RS485 interface have the following designation with respect to VP cable conductor colors:

| OXYFERM FDA ARC, OXYGOLD ARC | VP pin | VP8 cable |
|---|--------|-------------------------------|
| Power supply: +24 VDC (7 to 30 VDC). Max. power consumption 150 mW. | C | Coaxial core red transparent. |
| Power supply: Ground. | D | Coaxial shield red. |
| RS485 (A). | G | Yellow. |
| RS485 (B). | H | Brown. |

Table 3: RS485 interface pin layout

In an electromagnetically noisy environment, it is advisable to connect the VP cable shield to the ground. This significantly improves noise immunity and signal quality.

Examples of circuit arrangement


Fig. 6: Wiring diagram for the RS485 interface.


Fig. 7: Multi-drop bus wiring for the Modbus two-wire mode. Each sensor functions as a Modbus slave.

NOTE!

In the connection scheme shown above, each sensor must have the unique Modbus device address for proper communication.

The serial Modbus connection between the RS485 port of the master and the corresponding interfaces of the sensors has to be ensured according to the EIA/TIA RS485 standard. Only one sensor can communicate with the master at any time.

7. Configuration and monitoring of the sensor

Two options are available for the digital configuration and monitoring of a sensor:

1. Personal computer or notebook. The following additional equipment and software is also required:
 - HAMILTON USB-RS485 Modbus converter (P/N: 242411).
 - HAMILTON freeware 'ARC Sensor Configurator' available at <http://www.hamiltoncompany.com> (follow the instructions in the Configurator's user guide for the installation and operation).
 - Demo Cable (P/N: 355194). This cable includes a power adapter to supply the sensor with operation power, and a plug to connect the two RS485 conducting wires (yellow and brown) to USB-RS485. If you use a standard VP8 cable, you must supply the sensor with an external power source (pin C: 24 VDC; pin D: Ground).
2. HAMILTON ARC View Handheld Package (P/N: 242180). The handheld package represents an ideal desktop solution for the ARC sensor management. The ARC View Handheld included in the package is a compact mobile wireless device with long battery life and broad functionality. When using the Handheld as a mobile wireless device, each ARC sensor requires an ARC Wi Sensor Adapter (P/N: 242170).

The digital RS485 interface of the OXYFERM FDA ARC and OXYGOLD ARC sensor can be accessed by operators working within a three-tier hierarchy. Operator hierarchy levels and factory default passwords are shown in the table below.

| Operator status | Operator level | Password |
|-----------------|----------------|--------------|
| User | U | Not required |
| Administrator | A | 18111978 |
| Specialist | S | 16021966 |

Table 4: ARC operator levels

The Users can read basic data from sensor. Administrators can also calibrate sensors. Specialists can calibrate and configure sensors, and can see all data.

8. Preparation for the measurement

Prepare the sensor for measurements as follows:

1. Remove the protective cap from the VP head.
2. Make sure that the sensor is configured as required. If in doubt, test as described in the section entitled 'Electrical connection: RS485 Modbus interface'.
3. For OXYFERM FDA ARC, remove the protective watering cap from the membrane.
4. Unscrew the membrane cartridge and pre-rinse it with electrolyte. The OXYGOLD ARC sensors have no electrolyte in the membrane cartridge when they are new. Then fill it with 1.5 ml of electrolyte. Carefully screw the membrane cartridge tight again. Do not touch or

bend the glass body. The electrolyte is aggressive.

5. Mount the sensor by using the thread PG 13.5 / G1¹/₄". Do not use an upside down mounting.
6. Connect the sensor to power supply. The polarization will be started automatically and the corresponding warning message will appear. The polarization takes 2 hours for the OXYFERM FDA ARC and OXYGOLD G ARC and 24 hours for OXYGOLD B ARC sensors. This stabilization time is necessary to reach the specification of zero current. When the polarization is completed the sensor is ready for calibration.
7. To calibrate, follow the Section 11 'Testing and maintenance' in this manual.
8. If the sensor response in air is too small or sluggish, see the chapter about maintenance below.
9. Attach the sensor according to the section 'Electrical Connection' in the desired configuration (analog 4–20 mA interface, digital RS485 interface or both).
10. If the sensor was disconnected from power source only for a short time, it should be polarize 2 x the time of disconnection again before usage.

NOTE

Each time you disconnect the sensor, the polarization timer will be reset and a warning will appear for the standard duration of the polarization. If the disconnection was short and the zero current is stable, you can ignore this warning.

9. Removal of the sensor

Before removing the sensor, make sure that no process medium can be withdrawn thereby.

Separate the sensor thread and pull out the sensor. Do not turn the sensor at the socket head because otherwise the sensor cap loosens itself, and humidity reaches the interior of the sensor!

10. Sterilization, autoclaving, CIP purifications

OXYFERM FDA ARC and OXYGOLD ARC sensors are carefully designed to withstand the cleaning procedures normally used in biotechnological industries without the need for special precautions. Nevertheless, as is the case with all sensors, it is seen that frequent cleanings or sterilizations with steam or hot and aggressive solutions lead to shortened life span of the sensor. Although for a number of reasons (variations in composition of cleaning agent, temperature gradient and so on) it is not possible to give general life span figures.

The cable/connector contacts must be clean and dry before sensor is connected to the cable.

ATTENTION!

The sensor can communicate over digital RS485 interface up to a temperature of 130°C and up to 110°C (100°C for OXYGOLD B ARC) over analog 4-20 mA interface.

11. Testing and maintenance

The OXYFERM FDA ARC and OXYGOLD ARC sensors provide two kinds of sensor calibration: automatic standard calibration, and product calibration.

The automatic standard calibration is performed over RS485 interface using the ARC View Handheld or the ARC Sensor Configurator freeware (see the sections entitled 'Configuration and monitoring of the sensor' and 'Accessories'). For the product calibration use the ARC View Handheld, or the ARC Sensor Configurator freeware.

Automatic standard calibration

The OXYFERM FDA ARC and OXYGOLD ARC sensor is calibrated at two points: in air and in an oxygen-free environment. During calibration, the sensor examines the correctness and stability of the oxygen and temperature signals.

Calibration at point 1:

1. Set an appropriate operator level (Administrator or Specialist).
2. Immerse the sensor into an oxygen-free environment.
3. Let the system equilibrate. Guarantee stable conditions for at least three minutes. For greater measurement accuracy, make sure that temperature difference between calibration medium and process medium is minimal.
4. Execute the calibration at point 1.
5. If all conditions are met, the sensor confirms the calibration immediately. The calibration curve of the sensor is now defined by the recently calibrated point 1 and the prior calibration at point 2.

Calibration at point 2:

6. Set an appropriate operator level (Administrator or Specialist).
7. Leave the sensor for at least three minutes under stable conditions in ambient air or in oxygen saturated medium.
8. Execute the calibration at point 2.
9. If all conditions are met, the sensor confirms the calibration immediately. The calibration curve of the sensor is now defined by the recent calibration at both points: 1 and 2.

HINT: The concept behind OXYFERM FDA ARC and OXYGOLD ARC sensors enables calibration and configuration in the lab before use in the process control. Another calibration for the installation in the process setup is not required.

Product calibration.

The product calibration is an in-process calibration procedure in order to adjust the measurement to specific process conditions, or in case the sensor cannot be removed for the standard calibration.

The product calibration is an additional calibration procedure to a standard calibration. Product calibration adapts the standard calibration curve to the process conditions in force at the time of product calibration. If product calibration is activated, the calibration curve is calculated from the data of the last calibration at point 1 and from the data of the product calibration. In order to restore the original standard calibration curve, the product calibration can

be cancelled at any time. A new standard calibration cancels a product calibration as well.

Actions:

1. Set an appropriate operator level (Administrator or Specialist).
2. Perform an initial measurement while taking a sample from the process. The data of the initial measurement are stored in the sensor.
3. Perform a laboratory measurement of the sample at the same temperature as it was measured in the process.
4. Assign the laboratory value to the value of the initial measurement. This new DO value is accepted and instantaneous active, if the difference between initial measurement and laboratory values is not greater than 20%-sat units.

**NOTE!**

The product calibration is possible for DO values in the range of 2 to 55 %-vol (20 – 550 mbar).

Maintenance

Sensor check and maintenance for OXYFERM FDA ARC and OXYGOLD ARC sensor

For a periodical checking of the OXYGOLD ARC sensor please proceed as follows:

1. Hold a stabilized sensor in air and read measurement value. Hold the sensor in nitrogen flow and read measurement value (e.g. place nitrogen in a plastic bag). Measurement value must be under 2% of air value. Too high a current in nitrogen is due to an exhausted electrolyte, or a defective membrane. In such a case replace the membrane cartridge as well as the electrolyte solution.
2. If there is a sluggish response and/or too small a current in air, follow these steps:

| | | |
|---|-----|---------------------|
| Rinse connector with isopropanol. Dry thoroughly. | | |
| Specs ok? | yes | Reinsert in process |
| no | | |
| Clean membrane with water and a soft tissue | | |
| Specs ok? | yes | Reinsert in process |
| no | | |
| Exchange electrolyte | | |
| Specs ok? | yes | Reinsert in process |
| no | | |
| Replace membrane cartridge | | |
| Specs ok? | yes | Reinsert in process |
| no | | |
| Polish cathode (see warning below!) | | |
| Specs ok? | yes | Reinsert in process |
| no | | |
| Replace cathode | | |
| Specs ok? | yes | Reinsert in process |
| no | | |
| Contact HAMILTON dealer for support | | |

“Specs OK” means: current in air of stabilized sensor at 25°C is 200 ... 500 nA (OXYGOLD ARC) and 40 ... 80 nA (OXYFERM FDA ARC), $t_{98\%} < 60$ s at 25°C.


WARNING!

Do not polish the cathode of the OXYGOLD B ARC sensor!

Expanding the electrolyte and replacing the membrane cartridge

The OXYFERM FDA ARC and OXYGOLD ARC membrane cartridge is replaced as follows:

1. Put the sensor in an upright position and carefully unscrew the membrane cartridge.
2. Carefully clean the tip of the glass body with a soft tissue.
3. Check the small o-ring above the glass body. Replace if damaged.
4. Use the 1 ml plastic pipette of the membrane kit to pipette 1.5 ml of electrolyte (see below) into the new membrane cartridge.
5. Carefully screw the cartridge onto the sensor shaft. Any spillage of electrolyte should be rinsed away with water (see “Safety instructions”).

Polishing Cathode

The cathode and the front-end of the glass body need to be cleaned with the polishing strip included in the membrane kit. First take the dry polishing strip in one hand; then press the sensor lightly into the polishing strip and turn it 10-20 times. Afterwards, rinse the glass body under running water and dry it carefully.

ATTENTION

Mechanical polishing of the cathode should only be done with Hamilton polishing cloths. There is a polishing cloth included with the membrane kit. Do not touch the anode wire. Do not bend the glass body.

ATTENTION


Do not polish the cathode of an OXYGOLD B ARC sensor!

Replacing the Cathode

If the cathode is going to be changed, the membrane must also be changed to get the full performance of the sensor.

1. Put the sensor in an upright position and carefully unscrew the membrane cartridge.
2. Flush the replaceable cathode with deionised water then dry the metal parts. Do not touch anode and cathode.
3. Hold the cathode on the metal part (B) in front of the thread (A) and the sensor on the shaft (C) and pull apart. Just pull, do not turn!
4. Please check that all contacts are dry and clean.
5. Replace cathode for a new (p/n 237306 (blue) for OXYFERM FDA ARC; p/n 237427 (yellow) for OXYGOLD B ARC; p/n 237437 (black) for OXYGOLD G ARC). Rotate the cathode until the correct position is found and then push the cathode into the shaft.
6. Check the small O-ring (D) above the cathode. Replace if damaged.

7. Use the 1 ml plastic pipette of the membrane kit to pipette 1.5 ml of electrolyte into the new membrane cartridge. The pipette must NOT touch the membrane itself!
8. Screw the cartridge onto the sensor shaft. Any spillage of electrolyte should be rinsed away with water (see “Safety Instructions”).

Self-diagnosis functions

OXYFERM FDA ARC and OXYGOLD ARC sensors provide a self-diagnosis functionality to detect and identify the most common sensor malfunctions. Both interfaces can be used for warning and error messaging. The analog 4–20 mA interface can be configured according to the NAMUR recommendations to indicate an abnormal event. The RS485 interface provides a variety of the indications based on the fault code.

Use the ARC View Handheld and ARC Wi Sensor Adapter for monitoring of the sensor status and troubleshooting. These devices detect and highlight the sensor condition according to the fault code. The ARC View Handheld displays the text messages that are appropriate for the type of malfunction.

The following types of messages are provided by the self-diagnosis function:

- Warning (alarm):
 - Calibration recommended
 - Calibration upper/lower, out of range
 - Polarization stabilisation
- Error (failure):
 - DO sensor defect (reading failure)
 - Temperature sensor defect (reading failure)
 - Temperature upper/lower, out of range
 - Cathode impedance upper/lower range
 - DO partial pressure exceeds air pressure
 - Sensor quality low (calibration is not possible)
 - Internal communication failure

12. Disposal

The design of HAMILTON sensors optimally considers environmental compatibility. In accordance with the EC guideline 2002/96/EG HAMILTON sensors that are worn out or no longer required must be sent to a dedicated collection point for electrical and electronic devices, alternatively, must be sent to HAMILTON for disposal. Sensors must not be sent to an unsorted waste disposal point.

13. Parts and Accessories

ARC View Handheld with Dock.

ARC Wi Sensor Adapter for wireless communication.

USB-RS485 Modbus Converter.

Sensor Cable VP8.

| Order number | Description |
|---|--|
| For OXYFERM FDA ARC sensors | |
| 237123 | OXYFERM Membrane Kit |
| 237126 | CIP Membrane Kit |
| 237140 | FDA Membrane Kit |
| 237118 | OXYLYTE Electrolyte 50 ml |
| 237137 | Optional protective cap with chain |
| 237306 | Replacement Cathode OXYFERM |
| For OXYGOLD B ARC sensors | |
| 237135 | MEMBRANE KIT (3 cartridges + O-ring) |
| 237138 | OXYLYTE B Electrolyte 50 ml |
| 237427 | Replacement Cathode OXYGOLD B ARC |
| For OXYGOLD G ARC sensors | |
| 237135 | MEMBRANE KIT (3 cartridges + O-ring) |
| 237139 | OXYLYTE G Electrolyte 50 ml |
| 237437 | Replacement Cathode OXYGOLD G ARC |
| For all OXYFERM FDA ARC and OXYGOLD ARC sensors | |
| 355217 | Sensor cable VP8, 1m |
| 355218 | Sensor cable VP8, 3m |
| 355219 | Sensor cable VP8, 5m |
| 355220 | Sensor cable VP8, 10m |
| 355221 | Sensor cable VP8, 15m |
| 355222 | Sensor cable VP8, 20m |
| 242411 | USB-RS485 Modbus converter |
| 355194 | Demo Cable (1m, open end, with power pack) |
| 242180 | ARC View Handheld Package |
| 242170 | ARC Wi Sensor Adapter |

Table 5: Accessories.

14. Technical data

OXYFERM FDA ARC:

- ARC sensor with integrated electronics.
Functionality includes measurement and self-diagnosis.
- Electrode system: silver-platinum combination.
- Electrolyte: alkaline electrolyte solution.
- Membrane: OPTIFLOW.
- Shaft diameter: 12 mm standard; 25 mm XL-Type.
- Shaft length: several versions available.
- Operating temperature range:
 - Analog interface: 0 to 100°C.
 - Digital interface: 0 to 130°C.
- Storage temperature: -10 to 60 °C .
- Process pressure range: 0 to 4 bar.
- Detection limit: 10 ppb.
- Sensor current in air at 25 °C: 40 to 80 nA.
- Stabilizing time: max. 2 h.
- Mounting: PG 13.5 - thread standard; G 1¼" XL-Type.
- Response time $t_{98\%}$: max. 60 s at 25 °C, from air to nitrogen.
- Flow: min. 0.03 m/s.
- Flow response: < 5% at 25°C.
- Oxygen demand: ca. 20 ng/h in air at 25°C.
- Max. CO₂ partial pressure: 0.01 bar.
- Measured values can be configured by software according to:
 - DO: mbar; %-sat; %-vol; ug/l; mg/l; ppb; ppm.
 - Temperature: °C; °F; K.
- Materials: see certificate.
- Electrical connection: VP8 socket head.
- Operating voltage: 7 to 30 VDC max. 150 mW.
- Two scalable 4–20 mA current interfaces (two-wire) for DO and temperature signals.
- Digital RS485 Interface (two-wire):
 - Protocol: Modbus RTU; maximal 31 addresses.
 - Baud rate: 4800, 9600, 19200, 38400, 57600, 115200 bd.

OXYGOLD G ARC

- ARC sensor with integrated electronics.
Functionality includes measurement and self-diagnosis.
- Electrode system: silver-platinum combination.
- Electrolyte: alkaline electrolyte solution.
- Membrane: OPTIFLOW.
- Shaft diameter: 12 mm.
- Shaft length: several versions available.
- Operating temperature range:
 - Analog interface: 0 to 100°C.
 - Digital interface: 0 to 130°C.
- Storage temperature: -10 to 60 °C .
- Process pressure range: 0 to 12 bar.
- Detection limit: 1 ppb.
- Sensor current in air at 25 °C: 180 to 500 nA.
- Stabilizing time: max. 2 h.
- Mounting: PG 13.5 thread.
- Response time $t_{98\%}$: max. 60 s at 25 °C, from air to nitrogen.
- Flow: min. 0.1 m/s.
- Oxygen demand: ca. 100 ng/h in air at 25°C.
- Max. CO₂ partial pressure: 0.01 bar.
- Measured values can be configured by software according to:
 - DO: mbar; %-sat; %-vol; ug/l; mg/l; ppb; ppm.
 - Temperature: °C; °F; K.
- Materials: see certificate.
- Electrical connection: VP8 socket head.
- Operating voltage: 7 to 30 VDC max. 150 mW.
- Two scalable 4–20 mA current interfaces (two-wire) for DO and temperature signals.
- Digital RS485 Interface (two-wire):
 - Protocol: Modbus RTU; maximal 31 addresses.
 - Baud rate: 4800, 9600, 19200, 38400, 57600, 115200 bd.

OXYGOLD B ARC

- ARC sensor with integrated electronics.
Functionality includes measurement and self-diagnosis.
- Electrode system: silver-platinum combination.
- Electrolyte: acidic electrolyte solution.
- Membrane: OPTIFLOW.
- Shaft diameter: 12 mm.
- Shaft length: several versions available.
- Operating temperature range:
 - Analog interface: 0 to 100°C.
 - Digital interface: 0 to 130°C.
- Storage temperature: -10 to 60 °C .
- Process pressure range: 0 to 12 bar.
- Detection limit: 8 ppb.
- Sensor current in air at 25 °C: 180 to 500 nA.
- Stabilizing time: max. 24 h.
- Mounting: PG 13.5 - thread.
- Response time $t_{98\%}$: max. 60 s at 25 °C, from air to nitrogen.
- Flow: min. 0.1 m/s.
- Oxygen demand: ca. 100 ng/h in air at 25°C.
- Max. CO₂ partial pressure: 12 bar.
- Measured values can be configured by software according to:
 - DO: mbar; %-sat; %-vol; ug/l; mg/l; ppb; ppm.
 - Temperature: °C; °F; K.
- Materials: see certificate.
- Electrical connection: VP8 socket head.
- Operating voltage: 7 to 30 VDC max. 150 mW.
- Two scalable 4–20 mA current interfaces (two-wire) for DO and temperature signals.
- Digital RS485 Interface (two-wire):
 - Protocol: Modbus RTU; maximal 31 addresses.
 - Baud rate: 4800, 9600, 19200, 38400, 57600, 115200 bd.

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Copyright © 2011 HAMILTON Bonaduz AG, Bonaduz Schweiz. Alle Rechte vorbehalten. Die Reproduktion irgendeines Teils dieses Dokuments in jeder beliebigen Form ist ohne die ausdrückliche schriftliche Zustimmung der HAMILTON Bonaduz AG untersagt.