

Easy Arc Integration to new BioControllers

Technical Note

Standard M12-Pole Connectors fits best for new BioControllers

The M12 has become the industry standard in various applications. The wide selection of designs enables a cost efficient plug and play connection for different devices and sensors.

Hamilton has chosen the M12-8 Pole as a standard.




Figure 1: BioController scheme

A-Coded female connector

Pin	Function	
1	+ 4-20 mA	
2	- 4-20 mA	(5)
3	+ 4-20 mA	//(4) ⁽⁶⁾ (6)
4	- 4-20 mA	((3) (8) (7))
5	RS485 (A)	
6	RS485 (B)	
7	GND	
8	+ 24 VDC	Figure 2: M12 8 Pole A Code
Housing	Shield	female Connector

Arc Sensor cables for new BioControllers

Ref	Description	
10067842	3 m Sensor Cable VP8 / M12-8	
10067844	5 m Sensor Cable VP8 / M12-8	
10067846	10 m Sensor Cable VP8 / M12-8	

Arc Sensor Integration in to the BioController HMI System

We recommend to integrate these registers into the BioController HMI:

Register	Description	Read/Write
2090	Measured Value reading incl. sensor status (e.g. traffic light syml (DO, pH, COND, ORP)	bolic) Read
4736 (Warnings) 4800 (Errors)	Detailed Warning & Errors status	Read
4872	Quality indicator in % (DO, pH, COND, ORP)	Read
5158 (CP1) 5190 (CP2)	Calibration data Calibration Point 1 (DO, pH, COND, ORP) Status value: e.g. 4.01 pH Calibration point 2 (DO, pH) Status Value: e.g. 20.95 %-vol	Read
1280 (Ref) 1312 (SN) 1032 (FW)	Sensor ID Sensor ref number (Ref) Sensor serial number (SN) Sensor firmware (FW)	Read
4692 3146 3114 3530	Possible user input: (manual counter no input from the sensor) - Autoclaving counter (All parameter) - Atmospheric Pressure (DO only) - Salinity (DO only) - Cap Type (only for optical sensors only)	Write

These registers have to be set after each Arc-sensor has been powered up.

Register	Description	
4096	Device Address 1-32 (Default Address =1)	Write
4102	Baudrate (4800, 9600, 19200, 38400, 57600, 115200 bps)	Write
4288	Operator Level S	Write
2090	Measured Unit e.g. DO (define the specific parameter measured unit)	Write
4108	Parity	Write
4108	StopBits	Write

Table 1: Registers to be integrated into the BioController

Arc Standard Calibration:

Register	Description	
4288	Change Sensor to Operator Level S (Specialist)	Write
5162 (auto)	Start CP1 (Place Sensor in calibration media (wait at least 3 minutes for stable measurement and temperature condition)	Write
5194 (auto)	Start CP2 (Place Sensor in calibration media (wait at least 3 minutes for stable measurement and temperature condition)	Write
5158 (CP1) 5190 (CP2)	Read calibration Status for calibration points (CP1 or CP2)	Read

Arc Product Calibration:

Register	Description	
5322 (assign)	Assign offline a reference value	Write
5340 (init.)	Start Product Calibration with following values: Ox01 Perform initial measurement Ox02 Cancel an active product calibration	Write



Sensor and Bio Controller Management

Controller mode:

- Select Operator Level «S» (Exclusive access to the BioController)
- 2. Read and write rights via BioController (see table 1)
- 3. Write register after each sensor connection (see table 1)
- 4. Read out via wireless: everything that is read only («U»)
- 5. Sensor Calibration via BioController only
- 6. Product Calibration re-calibration via BioController only

Maintainance mode:

- 1. Select Operator Level «U» (unlock sensor for mobile or PC)
 - User has to select a button on the BioController HMI
 - Show user on HMI that read out via BioController is interrupted

2. New status:

- Read/Write out via BioController: not possible
- Read and write rights via ArcAir: everything

3. Switch to controller mode:

- User has to select button on BioController HMI
- When operator level is switched from level U to level «S» the following registers must be be overwritten by BioController:
 - Baud rate
 - Device address
 - Stop bits
 - Parity
 - Measurement unit




Figure 3: Controller mode




Figure 4: Maintainance mode

© 2018 Hamilton Bonaduz AG. All rights reserved.

REF 695246/00 — M 09/2018



Web: www.hamiltoncompany.com USA: 800-648-5950 Europe: +41-58-610-10-10

4970 Energy Way Reno, Nevada 89502 USA Tel: +1-775-858-3000 Fax: +1-775-856-7259 sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa

Via Crusch 8 CH-7402 Bonaduz, Switzerlan Tel: +41-58-610-10-10 Fax: +41-58-610-00-10 contact.pa.ch@hamilton.ch