

Transmitter H100 COND

User Manual



Latest Product Information:
www.hamiltoncompany.com



83529



HAMILTON 

Warranty

Defects occurring within 3 years from delivery date shall be remedied free of charge at our plant (carriage and insurance paid by sender).

Subject to change without notice.

Return of Products Under Warranty

Please contact our Service Team before returning a defective device. Request a Return Goods Authorization number (RGA). This number assures proper tracking of your device. Ship the cleaned device to the address you have been given.

If the device has been in contact with process fluids, it must be decontaminated/disinfected before shipment. In that case, please attach a corresponding certificate, for the health and safety of our service personnel.

Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

HAMILTON Bonaduz AG

Via Crusch 8

CH-7402 Bonaduz

Switzerland

Tel: +41 81 660 60 60

Fax: +41 81 660 60 70

Email: contact@hamilton.ch

Web: www.hamiltoncompany.com

HAMILTON Company

4970 Energy Way

Reno, Nevada 89502 USA

Tel: +1 775 858 3000

Toll Free: 800 648 5950

Fax: +1 775 856 7259

Email: sensors@hamiltoncompany.com

Web: www.hamiltoncompany.com

Hamilton Bonaduz AG Shanghai Office

Room 502, No.50 Boxia Road

Shanghai Pudong Software Park

201203 Shanghai

China

Phone: +86 21 61646567

Fax: +86 21 61063732

eMail: contact.china@hamilton.ch

Table of Contents

Safety Information	5
Intended Use.....	7
Registered Trademarks.....	7
CD-ROM	8
Safety Information	8
Quickstart Guides.....	8
Overview of Transmitter H100 COND.....	9
Assembly.....	10
Package Contents.....	10
Mounting Plan.....	11
Pipe Mounting, Panel Mounting.....	12
Installation and Connection.....	14
Installation Instructions	14
Terminal Assignments	14
Wiring of Conducell Sensors	17
Protective Wiring of Relay Outputs	18
User Interface and Display	20
Operation: Keypad	22
Safety Functions	23
Sensocheck, Sensoface Sensor Monitoring.....	23
GainCheck Device Self-Test.....	23
Automatic Device Self-Test	23
Hold Mode	24
Configuration	26
Menu Structure of Configuration.....	27
Overview of Configuration Steps	28
Output 1	30
Output 2	42
Temperature Compensation	48
Alarm Settings.....	50

Table of Contents

Limit Function.....	52
Controlling a Rinsing System.....	54
Connecting a Rinsing System.....	55
Parameters.....	56
Factory Settings of Parameters	56
Parameters – Individual Settings.....	58
Calibration	60
Calibration by Entry of Cell Constant.....	62
Calibration with Calibration Solution	64
Product Calibration	66
Temp Probe Adjustment	68
Measurement	68
Diagnostics Functions	69
Error Messages (Error Codes).....	71
Operating States.....	73
Sensoface	74
Appendix.....	77
Product Line and Accessories.....	77
Specifications.....	78
Calibration Solutions.....	83
Concentration Curves.....	86
Index	92
Passcodes.....	96

Safety information –

Be sure to read and observe the following instructions!

The device has been manufactured using state of the art technology and it complies with applicable safety regulations.

When operating the device, certain conditions may nevertheless lead to danger for the operator or damage to the device.

Caution!

Commissioning must be carried out by trained experts.

Whenever it is likely that protection has been impaired, the device shall be made inoperative and secured against unintended operation.

The protection is likely to be impaired if, for example:

- the device shows visible damage
- the device fails to perform the intended measurements
- after prolonged storage at temperatures above 70°C
- after severe transport stresses

Before recommissioning the device, a professional routine test in accordance with EN 61010-1 must be performed. This test should be carried out at the manufacturer's factory.

Caution!

Before commissioning, make sure that the transmitter may be connected to other equipment.

Intended Use

The Transmitter H100 COND is used for measurement of electrical conductivity and temperature in liquids. Fields of application are: biotechnology, chemical industry, environment, food processing, water/waste-water treatment.

The sturdy molded enclosure can be fixed into a control panel or mounted on a wall or at a post.

The protective hood provides additional protection against direct weather exposure and mechanical damage.

The device can be used with all 2- and 4-electrode sensors. It provides two current outputs (for transmission of measured value and temperature, for example), two contacts, and a universal power supply 24 ... 230 V AC/DC, AC: 45 ... 65 Hz.

Registered Trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

Sensocheck®

Sensoface®

GainCheck®

Provided Documentation



CD-ROM

Complete documentation:

- User manuals
- Safety instructions
- Quickstart guides



Safety Instructions

In official EU languages and others.

- EC Declarations of Conformity




Quickstart Guides

In German, English, French, Spanish.

- Installation and Commissioning
- Operation
- Menu structure
- Calibration
- Error messages and recommended actions

Overview of Transmitter H100 COND




Assembly

Package Contents

Check the shipment for transport damage and completeness.
The package should contain:

- Front unit
- Rear unit
- Bag containing small parts
- CD-ROM with documentation
- Specific test report
- Passcode sticker



- | | |
|---|--|
| 1 Jumper (2 x) | 9 Filler plug (3 x) |
| 2 Washer (1 x), for conduit mounting:
Place washer between enclosure and nut | 10 Hexagon nut (5 x) |
| 3 Cable tie (3 x) | 11 Sealing plug (2 x), for sealing in case of wall
mounting |
| 4 Hinge pin (1 x), insertable from either side | |
| 5 Enclosure screw (4 x) | |
| 6 Sealing insert (1 x) | |
| 7 Rubber reducer (1 x) | |
| 8 Cable gland (3 x) | |

Fig.: Assembling the enclosure

Mounting Plan





Fig.: Mounting plan (All dimensions in mm!)

Pipe Mounting, Panel Mounting



- 1 P/N 243084 protective hood (if required)
- 2 Hose clamp with worm gear drive to DIN 3017 (2 x)
- 3 Pipe-mount plate (1 x)
- 4 For vertical or horizontal posts or pipes
- 5 Self-tapping screw (4 x)

Fig.: P/N 243082 pipe-mount kit (All dimensions in mm!)





Fig.: P/N 243084 protective hood for wall and pipe mounting
(All dimensions in mm!)



- 1 Screw (4 x)
- 2 Gasket (1 x)
- 3 Control panel
- 4 Span piece (4 x)
- 5 Threaded sleeve (4 x)

Panel cut-out
138 x 138 mm (DIN 43700)

Fig.: P/N 243083 panel-mount kit (All dimensions in mm!)

Installation and Connection

Installation Instructions

Caution!

- Installation of the transmitter must be carried out by trained experts in accordance with this user manual and as per applicable local and national regulations.
- Be sure to observe the technical specifications and input ratings during installation.
- Be sure not to notch the conductor when stripping the insulation.
- Before connecting the device to the power supply, make sure that its voltage lies within the range 20.5 ... 253 V AC/DC.
- All parameters must be set by a system administrator prior to commissioning.

The terminals are suitable for single wires and flexible leads up to 2.5 mm² (AWG 14).

Terminal Assignments





Fig.: Transmitter H100 COND terminal assignments



- 1 Terminals for temperature probe and outer shield
- 2 Terminals for sensor
- 3 Terminals for power supply

Fig.: Information on installation, rear side of device

Wiring of Conducell Sensors

Wiring assignment for HAMILTON VP single coaxial cable, VP 6.0'

Transmitter H100 COND	Color coding in the cable	VP pin	Conducell 2 UP	Conducell 4 US Conducell 4 UxF
2	Coax core black/transparent	A	Pt ring 1	H Pot
1	Coax shield red	B	Pt ring 1	H Curr
3	Gray wire	C	Pt ring 2	L Pot
4	Blue wire	D	Pt ring 2	L Curr
D	White wire	E	Pt 1000	Pt 1000
E	Green wire	F	Pt 1000	Pt 1000
C	Outer shield green/yellow	Housing	Shielding of connector head ⁽¹⁾	Shielding of connector head ⁽¹⁾

Wiring assignment for HAMILTON VP double coaxial cable, VP 8.0 DC'


Transmitter H100 COND	Color coding in the cable	VP pin	Conducell 2 UP	Conducell 4 US Conducell 4 UxF
2	Coax core black/transparent	A	Pt ring 1	H Pot
1	Coax shield black	B	Pt ring 1	H Curr
3	Coax core red/transparent	C	Pt ring 2	L Pot
4	Coax shield red	D	Pt ring 2	L Curr
	White wire	E	Pt 1000	Pt 1000
	Green wire	F	Pt 1000	Pt 1000
D	Yellow wire	G	-	-
E	Brown wire	H	-	-
C	Outer shield green/yellow	Housing	Shielding of connector head ⁽¹⁾	Shielding of connector head ⁽¹⁾

(1) In the electrode the housing of the VP socket is not connected with any other potential-carrying component of the sensor.

Protective Wiring of Relay Outputs

Protective Wiring of Relay Contacts


Relay contacts are subjected to electrical erosion. Especially with inductive and capacitive loads, the service life of the contacts will be reduced. For suppression of sparks and arcing, components such as RC combinations, nonlinear resistors, series resistors and diodes should be used.



AC applications with inductive load

- 1 Load
- 2 RC combination, e.g. RIFA PMR 209
Typical RC combinations for 230 V AC:
Capacitor 0.1 μF / 630 V Resistor 100 ohms / 1 W
- 3 Contact

Typical Protective Wiring Measures



A: DC application with inductive load

B: AC/DC applications with capacitive load

C: Connection of incandescent lamps

A1 Inductive load

A2 Free-wheeling diode, e.g. 1N4007 (Observe polarity)

A3 Contact

B1 Capacitive load

B2 Resistor, e.g. $8\ \Omega$ / 1 W at 24 V / 0.3 A

B3 Contact

C1 Incandescent lamp, max 60 W / 230 V, 30 W / 115 V


C3 Contact

Warning!

Make sure that the maximum ratings of the relay contacts are not exceeded even during switching!


User Interface and Display

User Interface



- 1 Display
- 2 Mode indicators (no keys), from left to right:
 - Measuring mode
 - Calibration mode
 - Alarm
 - Cleaning contact
 - Configuration mode
- 3 Alarm LED
- 4 Keypad

Display



- | | | | |
|----|--|----|---------------------------|
| 1 | Passcode entry | 14 | Secondary display |
| 2 | Not in use | 15 | Manual temp specification |
| 3 | Temperature | 16 | Hold mode active |
| 4 | Current output | 17 | Waiting time running |
| 5 | Limit values | 18 | Sensor data |
| 6 | Alarm | 19 | Main display |
| 7 | Sensocheck | 20 | Sensoface |
| 8 | Calibration | | |
| 9 | Interval/response time | | |
| 10 | Cleaning contact | | |
| 11 | Measurement symbol | | |
| 12 | Press enter to proceed | | |
| 13 | Bar for identifying the device status,
above mode indicators, from left to right: | | |
| | - Measuring mode | | |
| | - Calibration mode | | |
| | - Alarm | | |
| | - Not in use | | |
| | - Configuration mode | | |

User Interface and Display

Operation: Keypad

cal	Start, exit calibration
conf	Start, exit configuration
▶	<ul style="list-style-type: none">• Select digit position (selected position blinks)• Menu navigation
▲	<ul style="list-style-type: none">• Edit digit• Menu navigation
enter	<ul style="list-style-type: none">• Calibration: Continue in program sequence• Configuration: Confirm entries, next configuration step• Measuring mode: Display output current

cal → enter	Cal Info, display of calibration data
conf → enter	Error Info: Display of last error message
▶ + ▲	Start GainCheck device self-test

Sensocheck, Sensoface Sensor Monitoring

Sensocheck continuously monitors the sensor and its wiring. Sensocheck can be switched off (Configuration, page 50).



Sensoface provides information on the conductivity sensor condition. Significant sensor polarization effects or an excessive cable capacitance are indicated.

GainCheck Device Self-Test

A display test is carried out, the software version is displayed, and the memory and measured-value transfer are checked.


Start GainCheck device self-test: ▶ + ▲

Automatic Device Self-Test

The automatic device self-test checks the memory and measured-value transfer. It runs automatically in the background at fixed intervals.

Safety Functions

Hold Mode

Display: 

The Hold mode is a safety state during configuration and calibration. Output current is frozen (Last) or set to a fixed value (Fix). Alarm and limit contacts are disabled.

If the calibration or configuration mode is exited, the device remains in the Hold mode for safety reasons. This prevents undesirable reactions of the connected peripherals due to incorrect configuration or calibration. The measured value and "HOLD" are displayed alternately. The device only returns to measuring mode after **enter** is pressed and 20 seconds have passed.

Configuration mode is also exited automatically 20 minutes (timeout) after the last keystroke. The device returns to measuring mode.

Timeout is not active during calibration.

Behavior of output signal:

- Last:** The output current is frozen at its last value.
Recommended for short configuration procedures.
The process should not change decisively during configuration.
Changes are not noticed with this setting!
- Fix:** The output current is set to a value that is noticeably different from the process value in order to signal the control system that the device is being worked at.

See Configuration page 40.

Alarm

Alarm delay is 10 seconds.




During an error message the alarm LED blinks.

Error messages can also be signaled by a 22 mA output current.

The alarm contact is activated by alarm or power failure, see also page 51.

Configuration

In the Configuration mode you set the device parameters.

<p>Activation</p>	<p>conf</p>	<p>Activate by pressing conf</p>
		<p>Enter passcode "1200" Edit parameter using ▶ and ▲, confirm/proceed using enter. (Exit by pressing conf, then enter.)</p>
<p>HOLD</p> <p>During configuration the device remains in the Hold mode.</p>	 <p>HOLD icon</p>	<p>The output current is frozen (at its last value or at a preset fixed value, depending on the configuration), limit and alarm contacts are inactive. Sensoface is off, "Configuration" mode indicator is on.</p>
<p>Input errors</p>		<p>The configuration parameters are checked during the input. In the case of an incorrect input "Err" is displayed for approx. 2 sec. The incorrect parameters cannot be stored. Input must be repeated.</p>
<p>Exit</p>	<p>conf</p> <p>enter</p>	<p>Exit by pressing conf. The measured value and Hold are displayed alternately, "enter" blinks. Press enter key to exit the Hold mode. The measured value is displayed. The output current remains frozen for another 20 sec (HOLD icon on, "hourglass" blinks).</p>

Menu Structure of Configuration

The configuration steps are assigned to different menu groups. Using the arrow keys, you can jump between the individual menu groups.

Each menu group contains menu items for setting the parameters. Pressing **enter** opens a menu item.

The values are edited using the arrow keys.

Pressing **enter** confirms/saves the settings.

Return to measurement: Press **conf**.

Select menu group	Menu group	Code	Display	Select menu item
	Output 1	o1.		enter
		Menu item 1		enter
		Menu item 2		enter
		:		enter
		Menu item ...		enter
	Output 2	o2.		
	Temperature compensation	tc.		
	Alarm settings	AL.		
	Relay	rL.		
	Rinsing system	Cn.		Previous menu group:

Configuration


Overview of Configuration Steps

Code	Menu	Selection
out1	Output 1	
o1.CELL	Select sensor	2-electrode, 4-electrode
o1.UnIT	Select process variable	μS, mS/cm, MΩ·cm, SAL, Conc
o1.CoNC	Select solution (Conc), see page 34	NaCl HCl NaOH H ₂ SO ₄ HNO ₃
	Codes:	-1- -2- -3- -4- -5-
o1.rNG	Select current range	0-20 mA / 4-20 mA
o1. 4mA	Enter current start	xxxx mS
o1.20mA	Enter current end	xxxx mS
o1.FtME	Time constant of output filter	xxxx SEC
o1.FAIL	22 mA signal in the case of error	ON / OFF
o1.HoLD	Signal behavior during HOLD	Last / Fix
o1.FIX	Enter fixed value	xxx.x mA
out2	Output 2	
o2.UnIT	Select temperature unit	°C / °F
o2. rTD	Select temperature probe	Pt100/Pt1000/NTC30 kΩ/ NTC8.55 kΩ
o2.rNG	Select current range	0-20 mA / 4-20 mA
o2. 4mA	Enter current start	xxx.x
o2.20mA	Enter current end	xxx.x
o2.FtME	Time constant of output filter	xxxx SEC
o2.FAIL	22 mA signal for temperature error	ON / OFF
o2.HoLD	Signal behavior during HOLD	Last / Fix
o2.FIX	Enter fixed value	xxx.x mA
tc.	Temperature compensation	
tc.	Select temp compensation	OFF/Lin/nLF/NaCl/HCl/NH ₃
tc. LIN	Lin: Enter temperature coefficient	xx.xx %/K

Code	Menu	Selection
ALrt	Alarm settings	
AL.SnSO	Select Sensocheck	ON / OFF
rLAY	Relay 1: Limit value	
L1.FCT	Select contact function	Lo / Hi
L1.tYP	Select contact response	N/O / N/C
L1.LEVL	Enter setpoint	xxxx
L1.HYS	Enter hysteresis	xxxx
L1.dLY	Enter delay	xxxx SEC
Cn	Rinsing probes	
Cn.InTV	Rinse interval	000.0 h
Cn.rins	Rinse duration	xxxx SEC
Cn.typ	Contact response	N/C / N/O


Configuration

Output 1 Selecting the sensor type




- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 **Output 1** menu group is displayed. All items of this menu group are indicated by the "o1." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 31). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.

Output 1:



4	
o1.CELL	Select sensor
o1.UnIT	Select process variable
o1.CoNC	Select solution (Conc)
o1.rNG	Select 0-20 / 4-20 mA
o1.4mA	Enter current start
o1.20mA	Enter current end
o1.FtME	Set output filter
o1.FAIL	22 mA for error
o1.HoLD	HOLD mode


5 **conf enter**

Code	Display	Action	Selection
01.		<p>Select evaluation method: 2-electrode sensor / 4-electrode sensor Select using ► key, press enter to proceed.</p>	<p>4-El (2-El / 4-El)</p>

Note: Characters represented in gray are blinking and can be edited.

Configuration

Output 1 Selecting the process variable



- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 **Output 1** menu group is displayed. All items of this menu group are indicated by the "o1." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 33). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.







Output 1:

3

4

o1.CELL	Select sensor
o1.UnIT	Select process variable
o1.CoNC	Select solution (Conc)
o1.rNG	Select 0-20 / 4-20 mA
o1.4mA	Enter current start
o1.20mA	Enter current end
o1.FtME	Set output filter
o1.FAIL	22 mA for error
o1.HoLD	HOLD mode


5 **conf enter**

Code	Display	Action	Choices
01.	     	<p>Select process variable:</p> <p>Select using ► key, press enter to proceed.</p> <p>Conductivity: 0.000 ... 9.999 µS/cm 00.00 ... 99.99 µS/cm 000.0 ... 999.9 µS/cm 0.000 ... 9.999 mS/cm 00.00 ... 99.99 mS/cm 000.0 ... 999.9 mS/cm 0.000 ... 9.999 S/m 00.00 ... 99.99 S/m</p> <p>Resistivity: 00.00 ... 99.99 MΩ·cm</p> <p>Salinity (SAL): 0.0 ... 45.0 ‰ (0 ... 35 °C)</p> <p>Concentration (Conc): 0.00 ... 9.99% by wt</p>	<p>000.0 mS (0.000 µS 00.00 µS 000.0 µS 0000 µS 0.000 mS 00.00 mS 000.0 mS 0.000 S/m 00.00 S/m 00.00 MΩ·cm 000.0 SAL 00.00 ‰)</p>

Configuration

Output 1

Concentration measurement: Select process solutions









- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 **Output 1** menu group is displayed. All items of this menu group are indicated by the "o1." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 35). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.

Output 1:

o1.CELL	Select sensor
o1.UnIT	Select process variable
o1.CoNC	Select solution (Conc)
o1.rNG	Select 0-20 / 4-20 mA
o1.4mA	Enter current start
o1.20mA	Enter current end
o1.FtME	Set output filter
o1.FAIL	22 mA for error
o1.HoLD	HOLD mode

5

conf enter

Code	Display	Action	Choices
o1.		Only with 00.00 % Conc can you select the process solution. Select using arrow key ▶	-01-SOL (-01-SOL -02-SOL -03-SOL -04-SOL -05-SOL)
		-01- NaCl (0.00 ... 9.99 % by wt) (0 ... 120 °C)	
		-02- HCl (0.00 ... 9.99 % by wt) (-20 ... 50 °C)	
		-03- NaOH (0.00 ... 9.99 % by wt) (0 ... 100 °C)	
		-04- H ₂ SO ₄ (0.00 ... 9.99 % by wt) (-17 ... 110 °C)	
		-05- HNO ₃ (0.00 ... 9.99 % by wt) (-20 ... 50 °C)	
		Press enter to proceed.	

Concentration Measurement


For the solutions listed above, the device can determine the substance concentration from the measured conductivity and temperature values in % by wt. The measurement error is made up of the sum of measurements errors during conductivity and temperature measurement and the accuracy of the concentration curves stored in the device, see page 86.

We recommend to calibrate the device together with the sensor. For exact temperature measurement, you should perform a temperature probe adjustment. For measuring processes with rapid temperature changes, use a separate temperature probe with fast response.

Configuration


Output 1

Output current range, current start, current end



- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 **Output 1** menu group is displayed. All items of this menu group are indicated by the "o1." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 37). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.

Output 1:



enter →




o1.CELL	Select sensor
o1.UnIT	Select process variable
o1.CoNC	Select solution (Conc)
o1.rNG	Select 0-20 / 4-20 mA
o1.4mA	Enter current start
o1.20mA	Enter current end
o1.FtME	Set output filter
o1.FAIL	22 mA for error
o1.HoLD	HOLD mode

4

enter

5



conf enter

Code	Display	Action	Choices
o1.		Set output current range Select using ▶ key, press enter to proceed.	4-20 mA (0 - 20 mA/ 4 - 20 mA)
		Current start Enter lower end of scale. Select using ▶ key, edit number using ▲ key, press enter to proceed.	000.0 mS (xxx.x mS)
		Current end Enter upper end of scale. Select using ▶ key, edit number using ▲ key, press enter to proceed.	100.0 mS (xxx.x mS)

Assignment of Measured Values: Current Start and Current End

Example 1: Range 0...200 mS/cm


Example 2: Range 100...200 mS/cm
Advantage: Higher resolution in the
range of interest



Configuration

Output 1

Time constant of output filter




- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 **Output 1** menu group is displayed. All items of this menu group are indicated by the "o1." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 39).
Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.

Output 1:

o1.CELL	Select sensor
o1.UnIT	Select process variable
o1.CoNC	Select solution (Conc)
o1.rNG	Select 0-20 / 4-20 mA
o1.4mA	Enter current start
o1.20mA	Enter current end
o1.FtME	Set output filter
o1.FAIL	22 mA for error
o1.HoLD	HOLD mode

5

conf enter

Code	Display	Action	Choices
o1.		Time constant of output filter Default setting: 0 s (inactive). To specify a time constant: Select using ► key, edit number using ▲ key, press enter to proceed.	0 sec 0 ... 120 sec

Time Constant of Output Filter (Attenuation)


To smoothen the current output, a low-pass filter with adjustable filter time constant can be switched on. When there is a jump at the input (100 %), the output level is at 63 % after the time constant has been reached.

The time constant can be set from 0 to 120 sec.

If the time constant is set to 0 sec, the current output follows the input.

Please note:


The filter only acts on the current output, not on the display or the limit value!



Configuration

Output 1

Output current during Error and HOLD



- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 **Output 1** menu group is displayed. All items of this menu group are indicated by the "o1." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 41). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.

Output 1:

3

enter

o1.CELL	Select sensor
o1.UnIT	Select process variable
o1.CoNC	Select solution (Conc)
o1.rNG	Select 0-20 / 4-20 mA
o1.4mA	Enter current start
o1.20mA	Enter current end
o1.FtME	Set output filter
o1.FAIL	22 mA for error
o1.HoLD	HOLD mode

4


enter

5

conf enter

Code	Display	Action	Choices
01.		22 mA signal for error message Select using ▶ key, press enter to proceed.	OFF (OFF/ON)
		Output signal during HOLD LAST: During HOLD the last measured value is maintained at the output FIX: During HOLD a value (to be entered) is maintained at the output Select using ▶ key, press enter to proceed.	LAST (LAST/FIX)
	 	Only with FIX selected: Enter current which is to flow at the output during HOLD Select position using ▶ key and edit number using ▲ key. Press enter to proceed.	21.0 mA (00.0 ... 21.0 mA)

Output Signal During HOLD:



Configuration

Output 2

Temperature unit and probe, output current

1 **conf**

2 1200

3 out. 1 MNU






Output 2:

3 out 2 MNU

4


o2.UnIT	Select °C/°F	enter
o2. rTD	Select temp probe	enter
o2.rNG	Select 0-20 / 4-20 mA	enter
o2. 4mA	Enter current start	
o2.20mA	Enter current end	
o2.FtME	Set output filter	
o2.FAIL	22 mA for temp error	
o2.HoLD	HOLD mode	

5 **conf** **enter**

Code	Display	Action	Choices
o2.		Specify temperature unit Select using ▶ key, press enter to proceed.	°C (°C/°F)
		Select temperature probe Select using ▶ key, press enter to proceed.	Pt1000 (Pt100, NTC30 kΩ, NTC8.55 kΩ)
		Select output current range Select using ▶ key, press enter to proceed.	4 - 20 mA (4 - 20 mA/ 0 - 20 mA)
		Current start: Enter lower end of scale. Select using ▶ key, edit number using ▲ key, press enter to proceed.	000.0 °C (xxx.x °C)
		Current start: Enter upper end of scale. Select using ▶ key, edit number using ▲ key, press enter to proceed.	100.0 °C (xxx.x °C)


Process Temperature: Current Start and Current End

Example 1: Range 0 ... 100 °C



Example 2: Range 50 ... 70 °C

Advantage: Higher resolution in the range of interest



Configuration

Output 2

Time constant of output filter


- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 Select **Output 2** menu group using arrow keys. All items of this menu group are indicated by the "o2." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 45). Confirm (and proceed) using **enter**.
Exit: Press **conf**, then **enter**.

Output 2:

o2.UnIT	Select °C/°F
o2. rTD	Select temp probe
o2.rNG	Select 0-20 / 4-20 mA
o2. 4mA	Enter current start
o2.20mA	Enter current end
o2.FtME	Set output filter
o2.FAIL	22 mA for temp error
o2.HoLD	HOLD mode

5

conf enter

Code	Display	Action	Choices
o2.		Time constant of output filter Default setting: 0 sec (inactive). To specify a time constant: Select using ▶ key, edit number using ▲ key, press enter to proceed.	0 sec (0 ... 120 sec)

Time Constant of Output Filter


To smoothen the current output, a low-pass filter with adjustable filter time constant can be switched on. When there is a jump at the input (100 %), the output level is at 63 % after the time constant has been reached.

The time constant can be set from 0 to 120 sec.

If the time constant is set to 0 sec, the current output follows the input.

Please note:

The filter only acts on the current output, not on the display!




Configuration


Output 2

Temperature error, output current during HOLD

1 **conf**



Output 2:



- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 Select **Output 2** menu group using arrow keys. All items of this menu group are indicated by the "o2." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 47). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.

enter





o2.UnIT	Select °C/°F
o2. rTD	Select temp probe
o2.rNG	Select 0-20 / 4-20 mA
o2. 4mA	Enter current start
o2.20mA	Enter current end
o2.FtME	Set output filter
o2.FAIL	22 mA for temp error
o2.HoLD	HOLD mode

enter


4

5

conf enter

Code	Display	Action	Choices
o2.		22 mA signal for error message Select using ▶ key, press enter to proceed.	OFF (OFF/ON)
		Output signal during HOLD LAST: During HOLD the last measured value is maintained at the output FIX: During HOLD a value (to be entered) is maintained at the output Select using ▶ key, press enter to proceed.	LAST (LAST/FIX)
	 	Only with FIX selected: Enter current which is to flow at the output during HOLD Select position with ▶ key and edit number with ▲ key. Press enter to proceed.	21.0 mA (00.0 ... 21.0 mA)






Output Signal During HOLD:










Configuration

Temperature Compensation

Temp compensation selection

- 1 **conf**
- 2 
- 3 

Temp compensation:

- 4 
- 5 **conf enter**


Code	Display	Action	Choices
tc.	     	<p>Select temp compensation</p> <p>OFF: Temperature compensation switched off Select using ► key, press enter to proceed.</p> <p>LIN: Linear temperature compensation with entry of temperature coefficient and reference temperature</p> <p>nLF: Temperature compensation for natural waters to EN 27888</p> <p>NaCl (nACL): Temperature compensation for ultrapure water with NaCl traces</p> <p>HCl (HCL): Temperature compensation for ultrapure water with HCl traces</p> <p>NH₃ (nH3): Temperature compensation for ultrapure water with NH₃ traces</p>	<p>OFF (OFF LIN nLF nACL HCL nH3)</p>
		<p>Only with linear temperature compensation (LIN) selected: Enter temperature coefficient. Select position using ► key and edit number using ▲ key. Press enter to proceed.</p>	<p>02.00%/K (XX.XX %/K)</p>

Configuration


Alarm Settings

1 **conf**


2




3



4




5



Alarm settings:

3




enter


4

AL.SnSO Select Sensocheck

5 **conf enter**

1 Press **conf** key.
2 Enter passcode **1200**.
3 Select **Alarm settings** menu group using arrow keys. All items of this menu group are indicated by the "AL." code.
4 Press **enter** to select menu, edit using arrow keys (see page 51). Confirm (and proceed) using **enter**.
5 Exit: Press **conf**, then **enter**.

Code	Display	Action	Choices
AL.		Select Sensocheck (continuous monitoring of sensor) Select using ► key, press enter to proceed.	OFF (ON/OFF)



Alarm Contact

The alarm contact is closed during normal operation (N/C). It opens in the case of alarm or power outage. As a result, a failure message is provided even in the case of line breakage (fail-safe behavior).

For contact ratings, see Specifications.


Error messages can also be signaled by a 22 mA output current (see page 40, 46, 71).

The operating behavior of the alarm contact is shown on page 73.

The **alarm delay** acts on the LED, the 22 mA signal and the alarm contact.

Configuration

Limit Function Relay




- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 Select **Limit function** menu group using arrow keys. All items of this menu group are indicated by the "L1." code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 53). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.

L1.FCT	Contact function	4 enter
L1.tYP	Contact response	
L1.LEVL	Enter setpoint	
L1.HYS	Enter hysteresis	
L1.dLY	Delay	


5 **conf enter**

Code	Display	Action	Choices
L1.		Contact function (see below for function principle) Select using ▶ key, press enter to proceed.	Lo (Lo/Hi)
		Contact response N/C: normally closed contact N/O: normally open contact Select using ▶ key, press enter to proceed.	N/O (N/O N/C)
		Setpoint Select using ▶ key, edit number using ▲ key, press enter to proceed.	000.0 mS (xxx.x mS)
		Hysteresis Select using ▶ key, edit number using ▲ key, press enter to proceed.	001.0 mS (xxx.x mS)
		Delay The contact is activated with delay (deactivated without delay) Select using ▶ key, edit number using ▲ key, press enter to proceed.	0010 sec (0 ... 9999 sec)


Limit Lo



Limit Hi






Controlling a Rinsing System “Clean” contact



- 1 Press **conf** key.
- 2 Enter passcode **1200**.
- 3 Select **Limit function** menu group using arrow keys. All items of this menu group are indicated by the “Cn.” code.
- 4 Press **enter** to select menu, edit using arrow keys (see page 54). Confirm (and proceed) using **enter**.
- 5 Exit: Press **conf**, then **enter**.


Cn.InTV	Rinsing interval	4 enter enter enter
Cn.rins	Rinse duration	
Cn.typ	Contact response	

5 conf enter

Code	Display	Action	Choices
Cn.		Rinsing interval Select using ▶ key, enter number using ▲ , press enter to proceed.	000.0 h (xxx.x h)
		Rinse duration Select using ▶ key, enter number using ▲ , press enter to proceed.	0060 sec (xxxx sec)
		Contact response N/C: normally closed contact N/O: normally open contact Select using ▶ , press enter to proceed.	N/C (N/O)

Connecting a Rinsing System

The “Clean” contact can be used to connect a simple rinsing system. Rinse duration and rinsing interval are defined during configuration.



Parameters

Factory Settings of Parameters

Activation:

Simultaneously press **conf** + right arrow key and enter passcode "4321".

The lower display line reads "Clear". To prevent accidental resetting, "NO" is set as default (blinking in the main display).

Press one of the arrow keys to select "YES" and confirm by pressing **enter**.

Caution!

Your data (also calibration data) will be overwritten by the factory settings!

Code	Parameters	Factory setting
o1.CELL	Select sensor	4-EI
o1.UnIT	Process variable	000.0 mS
o1.CoNC	Conc solution	-01-
o1.rNG	0/4-20 mA	4-20 mA
o1.4mA	Current start	000.0 mS
o1.20mA	Current end	100.0 mS
o1.FtME	Filter time	0 s
o1.FAIL	22mA signal	OFF
o1.HoLD	HOLD response	Last
o1.FIX	Fix current	021.0 mA
o2.UnIT	Unit °C / °F	°C
o2.rTD	Temp probe	Pt 1000
o2.rNG	0/4-20 mA	4-20 mA
o2.4mA	Current start	000.0 °C
o2.20mA	Current end	100.0 °C
o2.FtME	Filter time	0 s
o2.FAIL	22mA signal	OFF
o2.HoLD	HOLD response	Last
o2.FIX	Fix current	021.0 mA

Code	Parameters	Factory setting
tc.	Temp compensation	OFF
tc. LIN	Temp coefficient	02.00%/K
AL.SnSO	Sensocheck	OFF
L1.FCT	Contact function	Lo
L1.tYP	Contact response	N/O
L1.LEVL	Setpoint	000.0 mS
L1.HYS	Hysteresis	001.0 mS
L1.dLY	Delay	0010 sec
Cn.InTV	Rinsing interval	000.0 h
Cn.rins	Rinse duration	0060 sec
Cn.typ	Contact type	N/C

Please note:

Fill in your configuration data on the following pages.

Please note:

The cell constant is factory set to 1.0000 cm⁻¹.

Parameters






Parameters – Individual Settings

Code	Parameter	Setting
o1.CELL	Sensor	
o1.UnIT	Process variable	
o1.CoNC	Solution (Conc)	
o1.rNG	0/4-20 mA	
o1.4mA	Current start	
o1.20mA	Current end	
o1.FtME	Filter time	
o1.FAIL	22mA signal	
o1.HoLD	HOLD response	
o1.FIX	Fix current	
o2.UnIT	Unit °C / °F	
o2.rTD	Temp probe	
o2.rNG	0/4-20 mA	
o2.4mA	Current start	
o2.20mA	Current end	

Code	Parameter	Setting
o2.FtME	Filter time	
o2.FAIL	22mA signal	
o2.HoLD	HOLD response	
o2.FIX	Fix current	
tc.	Temp compensation	
tc. LIN	Temp coefficient	
AL.SnSO	Sensocheck	
L1.FCT	Contact function	
L1.tYP	Contact response	
L1.LEVL	Setpoint	
L1.HYS	Hysteresis	
L1.dLY	Delay	
Cn.InTV	Rinsing interval	
Cn.rins	Rinse duration	
Cn.typ	Contact type	

Calibration

Calibration adjusts the device to the sensor.

Activation	cal	Activate by pressing cal
		<p>Enter passcode:</p> <ul style="list-style-type: none"> • Entry of cell constant 1100 • With calibration solution 0110 • Product calibration 1105 • Temp probe adjustment 1015 <p>Select using ▲ key. Edit parameter using ▶ . Press enter to proceed. (Exit by pressing cal, then enter.)</p>
<p>HOLD</p> <p>During calibration the device remains in the Hold mode.</p>	   HOLD icon	<p>Output current is frozen (last value or preset fixed value, depending on configuration), limit and alarm contacts are inactive. Sensoface is off, "Calibration" mode indicator is on.</p>
<p>Input errors</p>		<p>The calibration parameters are checked during the input. In the case of an incorrect input "Err" is displayed for approx. 2 sec. The incorrect parameters cannot be stored. Input must be repeated.</p>
<p>Exit</p>	<p>enter</p> <p>enter</p>	<p>Exit by pressing enter (abort using cal).</p> <p>The measured value and Hold are displayed alternately, "enter" blinks. Sensoface is active. Press enter to exit the Hold mode. The measured value is displayed. The output current remains frozen for another 20 sec (HOLD icon on, "hourglass" blinks).</p>

Information on Calibration

Calibration adapts the device to the conductivity sensor.

Calibration can be performed by:

- Input of cell constant (e.g. for ultrapure-water sensors)
- Determining the cell constant with a known calibration solution (conductivity standard)
- Product calibration (calibration by comparison)
- Temperature probe adjustment







Please note:


- All calibration procedures must be performed by trained personnel.
- Incorrectly set parameters may go unnoticed, but change the measuring properties.

Calibration

Calibration by Entry of Cell Constant

Input of cell constant with simultaneous display of uncorrected conductivity value and temperature






Display	Action	Remark
	Press cal key, enter code 1100. Select using ▶ key, edit number using ▲ key, press enter to proceed.	Device is in the Hold mode. If an invalid code is entered, the device returns to measuring mode.
	Ready for calibration	Display (2 sec)
   	Enter cell constant of connected sensor: Select using ▶ , enter number using ▲ A change in the cell constant also changes the conductivity value. Press enter to confirm cell constant.	The lower display shows the measured conductivity value. (When there has not been an entry for 6 sec, the lower display alternately shows the conductivity and temperature value.)


Display	Action	Remark
 <p>The image shows a digital display with a smiley face icon in the top left. The main display shows '10.03 mS'. Below this, there is a 'hold' icon (a triangle with 'hold' written below it) and the temperature '26.3 °C'. To the right of the temperature, there is a small icon of a keypad.</p>	<p>The device now displays the conductivity and temperature.</p> <p>The measured value is shown in the main display alternately with "Hold", "enter" blinks. Exit calibration by pressing enter.</p>	<p>After end of calibration, the outputs remain in Hold mode for approx. 20 sec.</p>

Calibration

Calibration with Calibration Solution

Input of temperature-corrected value of calibration solution (calibration standard) with simultaneous display of cell constant

Display	Action	Remark
	Press cal key, enter code 0110. Select using ▶ key, edit number using ▲ key, press enter to proceed.	Device is in the Hold mode. If an invalid code is entered, the device returns to measuring mode.
	Ready for calibration Remove and clean sensor	Display (2 sec)
	Immerse sensor in calibration solution. Determine the temperature-corrected conductivity value of the calibration solution from the corresponding table (see page 83).	When there has not been an entry for 6 sec, the lower display alternately shows the cell constant and temperature value.
 	Enter value of calibration solution. Select using ▶ key, edit number using ▲ key. Press enter to confirm the calibration.	The cell constant and temperature are alternately displayed in the lower display during the input.
	The determined cell constant is displayed. Press enter to confirm.	

Display	Action	Remark
	<p>The device now displays the conductivity and temperature.</p>	
	<p>Clean sensor and re-place it in the process. The measured value is shown in the main display alternately with "Hold". "enter" blinks. Exit calibration by pressing enter.</p>	<p>After end of calibration, the outputs remain in Hold mode for approx. 20 sec.</p>

Please note:

- Be sure to use known calibration solutions with the respective temperature-corrected conductivity values (see "Calibration Solutions" page 83 et seq.).
- Make sure that the temperature does not change during the calibration procedure.





Calibration


Product Calibration

Calibration by comparison




For product calibration the measured variable is used as configured: conductivity ($\mu\text{S}/\text{cm}$, mS/cm , S/m), resistivity ($\text{M}\Omega\text{-cm}$). During product calibration the sensor remains in the process. The measurement is only interrupted briefly. Calibration is without TC correction.

Procedure: The currently measured value is stored in the device for comparison. A sample is measured using a portable meter. The sample value is then entered in the device. The new cell constant is calculated from these two values.


Display	Action	Remark
 The display shows the code '1105' in large digits. Above the first digit is a small lock icon. Below the digits are two small icons: a triangle with a lightning bolt and a square with a right-pointing arrow.	Press cal key, enter code 1105. Press ▶ key to select position, enter number using ▲ key, press enter to confirm.	If an invalid code is entered, the device returns to measuring mode.
 The display shows 'CAL PRD' in large letters. Above the 'A' is a small icon of a cat's head. Below the text is a horizontal line.		Display (approx. 2 sec)
 The display shows '1390 mS' in large digits. Above the '0' is a small icon of a cat's head. Below the digits is the word 'STORE' in smaller letters. At the bottom are two small icons: a triangle with a lightning bolt and a square with a right-pointing arrow.	Save currently measured value. Press enter to proceed.	Perform reference measurement.
 The display shows '1285 mS' in large digits. Above the '5' is a small icon of a cat's head. Below the digits is the word 'CALC' in smaller letters. At the bottom are two small icons: a triangle with a lightning bolt and a square with a right-pointing arrow.	Enter reference value. The new cell constant is calculated.	

Display	Action	Remark
	<p>The determined cell constant is displayed. Press enter to confirm.</p>	<p>New calibration: Press cal.</p>
	<p>The new value is shown in the main display alternately with "Hold", "enter" blinks. Exit by pressing enter.</p>	<p>After end of calibration, the outputs remain in Hold mode for approx. 20 sec.</p>





Temp Probe Adjustment

Display	Action	Remark
	Select calibration Press cal key, enter code 1015. Press ▶ key to select position, enter number using ▲ key, press enter to confirm.	Wrong settings change the measurement properties! If an invalid code is entered, the device returns to measuring mode.
	Ready for calibration Measure the temperature of the process medium using an external thermometer.	Device is in Hold mode. Display approx. 2 sec
	Enter measured temperature value. Select using ▶ key, edit number using ▲ key, press enter to proceed. Exit adjustment by pressing enter . HOLD will be deactivated after 20 sec.	Default: Value of secondary display.

Measurement





Display	Action
	In the measuring mode the main display shows the configured process variable (conductivity, resistivity, or SAL) and the lower display the temperature. During calibration you can return to measuring mode by pressing the cal key, during configuration by pressing conf (waiting time for signal stabilization approx. 20 sec).

Diagnostics Functions

Display	Action
 The display shows a main reading of 13.2 mA and a secondary reading of 12.5 mA. A small icon of a circle with a dot is in the top left corner.	Display of output currents Press enter while in measuring mode. The current at output 1 is shown in the main display, the current at output 2 in the secondary display. After 5 sec the device returns to measuring mode.
 The display shows 10.0000 in the main area and CELL in the secondary area. A smiley face icon is in the top left, and a small graph icon with the word SLOPE is in the bottom left.	Display of calibration data (Cal Info) Press cal while in measuring mode and confirm code 0000. The current cell constant is shown in the main display. After 20 sec the device returns to measuring mode (immediate return at pressing enter).
 The display shows 1002 KR in the main area and 38.2 OC in the secondary area. A small icon of a square with a diagonal line is in the bottom right.	Sensor monitor for validation of sensor and complete signal processing. Press conf while in measuring mode and enter code 2222. The measured resistance is shown in the main display, the measuring temperature in the lower display. Press enter to return to measurement.
 The display shows the word LAST in the main area and Error in the secondary area. A smiley face icon is in the top left, and a small icon of a square with a diagonal line is in the bottom left.	Display of last error message (Error Info) Press conf while in measuring mode and confirm code 0000. The last error message is displayed for approx. 20 sec. After that the message will be deleted (immediate return to measurement at pressing enter).

Diagnosics Functions















These functions are used for testing the connected peripherals.

Display	Action
	Specify current for output 1 Press conf while in measuring mode, enter code 5555. The current indicated in the main display for output 1 can be edited.
	Select using ▶ key, edit number using ▲ key. Press enter to confirm entry. The entered value will be shown in the secondary display. The device is in Hold mode. Press conf , then enter to return to measurement (Hold remains active for another 20 sec).
	Specify current for output 2 Press conf while in measuring mode, enter code 5556. The current indicated in the main display for output 2 can be edited.
	Select using ▶ key, edit number using ▲ key. Press enter to confirm entry. The entered value will be shown in the secondary display. The device is in Hold mode. Press conf , then enter to return to measurement (Hold remains active for another 20 sec).































Error Messages (Error Codes)

Error	Display	Problem Possible causes	Alarm contact	Red LED	Out 1 (22 mA)	Out 2 (22 mA)
ERR 01	Measured value blinks	Sensor <ul style="list-style-type: none"> • Wrong cell constant • Measuring range violation • SAL > 45 ‰ • Sensor connection or cable defective 	x	x	x	
ERR 02	Measured value blinks	Unsuitable sensor Conductance range > 3500 mS	x	x	x	
ERR 98	“Conf” blinks	System error Configuration or calibration data defective; completely reconfigure the device using the factory settings. Then calibrate. Memory error in device program	x	x	x	x
ERR 99	“FAIL” blinks	Factory settings EEPROM or RAM defective This error message only occurs in the case of a total defect. The device must be repaired and recalibrated at the factory.	x	x	x	x

Error Messages (Error Codes)

Error	Icon (blinks)	Problem Possible causes	Alarm contact	Red LED	Out 1 (22 mA)	Out 2 (22 mA)
ERR 03		Temperature probe Open or short circuit Temperature range exceeded	x	x	x	x
ERR 11		Current output 1 Current below 0 (3.8) mA	x	x	x	
ERR 12		Current output 1 Current above 20.5 mA	x	x	x	
ERR 13		Current output 1 Current span too small / too large	x	x	x	
ERR 21	 	Current output 2 Current below 0 (3.8) mA	x	x		x
ERR 22	 	Current output 2 Current above 20.5 mA	x	x		x
ERR 23	 	Current output 2 Current span too small / too large	x	x		x
ERR 33	 	Sensocheck: Wrong or defective sensor / polarization effects at the sensor / cable too long or defective / plug defective	x	x	x	
	 	Temperature outside conversion tables (TC, conc, SAL)				

Operating States

Operating status	Out 1	Out 2	Relay 1 limit value	Alarm contact	Cleaning contact	Timeout
Measure						
Cal Info (cal) 0000						20 s
Error Info (conf) 0000						20 s
Calibration (cal) 1100						
Temp adjustment (cal) 1015						
Product calibration (cal) 1105						
Configuration (conf) 1200						20 min
Sensor monitor (conf) 2222						20 min
Current source 1 (conf) 5555						20 min
Current source 2 (conf) 5556						20 min
Rinsing function						

 active

 as configured (Last/Fix or Last/Off)

Sensoface

The smiley in the display (Sensoface) provides information about the sensor condition (defects, maintenance required, cable capacitance too high). It alerts to significant sensor polarization or excessive cable capacitance e.g. caused by an unsuitable cable or a cable that is too long. The permitted calibration ranges and the conditions for a friendly, neutral, or sad Sensoface are summarized in the following chart. Additional icons refer to the error cause.

Sensocheck

Continuously monitors the sensor and its wiring. Sensocheck can be switched off. Critical values make the Sensoface “sad” and the corresponding icon blinks:







The Sensocheck message is also output as error message Err 33. The alarm contact is active, the red LED is lit, output current 1 is set to 22 mA (when configured correspondingly). Sensocheck can be switched off during configuration (then Sensoface is also disabled).

Exception: After a calibration a smiley is always displayed for confirmation.

Notice

The worsening of a Sensoface criterion leads to the devaluation of the Sensoface indicator (Smiley becomes “sad”). An improvement of the Sensoface indicator can only take place after calibration or removal of the sensor defect.

Display	Problem	Status
	Sensor defect	 Wrong or defective sensor Significant polarization of sensor Excessive cable capacitance (see also Err 33, Error Messages on page 72).
	Temperature error	 Temperature outside range for TC, conc, SAL

Please note:

When very fast response times (t_{90}) are required, e.g. when detecting separation layers, Sensoscheck should be switched off (see “Specifications” page 78).

Product Line and Accessories

Devices

Transmitter H100 COND

Part No.

243080-02

Mounting Accessories

Pipe-mount kit

243082

Panel-mount kit

243083

Protective hood

243084

Specifications

Conductivity input	Input for 2-electrode/4-electrode sensors	
Effective range	Conductivity	0.2 $\mu\text{S} \cdot \text{cm}$... 1000 $\text{mS} \cdot \text{cm}$
Measuring ranges	Conductivity	0.000 ... 9.999 $\mu\text{S}/\text{cm}$
		00.00 ... 99.99 $\mu\text{S}/\text{cm}$
		000.0 ... 999.9 $\mu\text{S}/\text{cm}$
		0000 ... 9999 $\mu\text{S}/\text{cm}$
		0.000 ... 9.999 mS/cm
		00.00 ... 99.99 mS/cm
		000.0 ... 999.9 mS/cm
		0,000 ... 9.999 S/m
		00.00 ... 99.99 S/m
		000.0 ... 999.9 S/m
	Resistivity	00.00 ... 99.99 $\text{M}\Omega \cdot \text{cm}$
	Concentration	0.00 ... 9.99 % by wt
	Salinity	0.0 ... 45 ‰ (0 ... 35 °C)
Response time (T_{90})	< 1 s (Sensocheck off)	
	< 3 s (Sensocheck on)	
Meas. error ^{1,2,3)}	< 1 % meas. val. + 0.4 $\mu\text{S} \cdot \text{cm}$	
Concentration determination		
Operating modes *	-01-	NaCl 0.00...9.99 % by wt (0...60 °C)
	-02-	HCl 0.00...9.99 % by wt (-20...50 °C)
	-03-	NaOH 0.00...9.99 % by wt (0...100 °C)
	-04-	H ₂ SO ₄ 0.00...9.99 % by wt (-17...110 °C)
	-05-	HNO ₃ 0.00...9.99 % by wt (-20...50 °C)
	See graphs in the Appendix, page 86	
Sensor standardization		
Operating modes	<ul style="list-style-type: none"> • Input of cell constant with simultaneous display of conductivity and temperature • Input of conductivity of calibration solution with simultaneous display of cell constant and temperature • Product calibration • Temperature probe adjustment 	
Permissible cell constant	00.0050 ... 19.9999 cm^{-1}	

Sensor monitoring

Sensocheck

Polarization detection and monitoring of cable capacitance

Sensoface

Provides information on the sensor condition (Sensocheck)

Sensor monitor

Direct display of measured values from sensor for validation (resistance/temperature)

Temperature input *

Pt100/Pt1000/NTC 30 k Ω /NTC 8.55 k Ω (Betatherm)
2-wire connection, adjustable

Measuring range

Pt 100/Pt 1000	-20.0 ... +200.0 °C (-4...+392 °F)
NTC 30 k Ω	-20.0 ... +150.0 °C (-4...+302 °F)
NTC 8.55 k Ω	-10.0 ... +130.0 °C (+14...+266 °F)

Resolution

0.1 °C / 0.1 °F

Meas. error^{1,2,3)}

< 0.5 K (< 1K for Pt100; < 1K for NTC > 100°C)

Temperature compensation

(reference temp 25°C)

(OFF)	Without
(Lin)	Linear characteristic 00.00 ... 19.99 %/K
(NLF)	Natural waters to EN 27888
(nACL)	Ultrapure water with NaCl traces (0...120°C)
(HCL)	Ultrapure water with HCl traces (0...120°C)
(nH3)	Ultrapure water with NH ₃ traces (0...120°C)

Output 1

0/4 ... 20 mA, max. 10 V, floating
(galvanically connected to output 2)

Process variable*

Conductivity, resistivity, concentration, salinity

Overrange*

22 mA in the case of error messages

Output filter*

Low-pass, filter time constant 0 ... 120 s

Measurement error¹⁾

< 0.3% current value + 0.05 mA

Start/end of scale

As desired within range

Minimum span

5 % of selected range

Specifications

Output 2	0/4 ... 20 mA, max. 10 V, floating (galvanically connected to output 1)
Process variable	Temperature
Overrange *	22 mA in case of temp error messages
Output filter *	Low-pass, filter time constant 0 ... 120 s
Measurement error ¹⁾	< 0.3% current value + 0.05 mA
Start/end of scale *	-20 ... 300 °C / -4 ... 572 °F
Permissible span	20 ... 320 K / 36 ... 576 °F
Alarm contact	Relay contact, floating
Contact ratings	AC < 250 V / < 3 A / < 750 VA DC < 30 V / < 3 A / < 90 W
Contact response	N/C (fail-safe type)
Alarm delay	10 s
Limit values	Output via relay contact
Contact ratings	AC < 250 V / < 3 A / < 750 VA DC < 30 V / < 3 A / < 90 W
Contact response *	N/O or N/C
Delay *	0000 ... 9999 s
Setpoints *	As desired within range
Hysteresis *	0 ... 50 % full scale
Rinsing function	Relay contact, floating, for controlling a simple rinsing system
Contact ratings	AC < 250 V / < 3 A / < 750 VA DC < 30 V / < 3 A / < 90 W
Contact response	N/C or N/O
Rinse interval	000.0 ... 999.9 h (000.0 h = cleaning function switched off)
Rinse duration	0000 ... 1999 s
Display	LC display, 7-segment with icons
Main display	Character height 17 mm, unit symbols 10 mm
Secondary display	Character height 10 mm, unit symbols 7 mm
Sensoface	3 status indicators (friendly, neutral, sad face)
Mode indication	4 mode indicators "meas", "cal", "alarm", "config" Further icons for configuration and messages
Alarm indication	Red LED in case of alarm

Keypad	5 keys: [cal] [conf] [▶] [▲] [enter]
Service functions	
Current source	Current specifiable for output 1 and 2 (00.00 ... 22.00 mA)
Device self-test	Automatic memory test (RAM, FLASH, EEPROM)
Display test	Display of all segments
Last Error	Display of last error occurred
Sensor monitor	Display of direct sensor signal (resistance/temperature)
Data retention	Parameters and calibration data > 10 years (EEPROM)
Protection against electric shock	Protective separation of all extra-low-voltage circuits against mains by double insulation to EN 61010-1
Power supply	24 (-15%)...230 V AC/DC (+10%); appr. 5 VA, 2.5 W AC: 45 ... 65 Hz Overvoltage category II, protection class II
Nominal operating conditions	
Ambient temperature	-20 ... +55 °C
Transport/Storage temp	-20 ... +70 °C
Relative humidity	80 % at temperatures up to 55 °C, maximum operating height 2000 m
Power supply	24 (-15%) ... 230 V AC/DC (+10%)
Frequency for AC	45 ... 65 Hz
EMC	EN 61326-1, EN 61326-2-3
Emitted interference	Class B (residential area) Class A for mains > 60 V DC
Immunity to interference	Industry

Specifications

Housing	Molded enclosure made of PBT (polybutylene terephthalate)
Color	Bluish gray RAL 7031
Mounting	<ul style="list-style-type: none">• Wall mounting• Pipe mounting: Ø 40 ... 60 mm □ 30 ... 45 mm• Panel mounting, cutout to DIN 43 700 Sealed against panel
Dimensions	H 144 mm, W 144 mm, D 105 mm
Ingress protection	IP 65 / NEMA 4X
Cable glands	3 knockouts for cable glands M20x1.5 2 knockouts for NPT 1/2" or rigid metallic conduit
Weight	Approx.1 kg

* User-defined

1) To IEC 746 Part 1, at nominal operating conditions

2) ± 1 count

3) Plus sensor error

Calibration Solutions

HAMILTON conductivity standards



DANAK Calibration Ltd.
Newcastle Business Park
Newcastle upon Tyne
Tel: +44 (0) 191 261 1111
www.danak.com

Certificate no. 02662
Page 1 of 2

Part No: **CON-113-0**
 Part No: **CON-113-0**
 Part No: **CON-113-0**

Calibration certificate

Electrolytic conductivity

Client	Macheca Biotech Ltd
Address	19a Church St, 1142 Hendon, London
Telephone/Fax	+41 61 848 6986
Contact person	Dr. Michael Maynard
Order reference	2008-08-13
Identification	Conductivity standard 1,0 µS/cm
Batch	R/N 2.00071, W/L 1345523
Date of calibration	2008-08-13

Result: Conductivity standard 1,0 µS/cm, P/N 23887A, W/L 1345523, Sample 3
 Laboratory conductivity standard (± 2.0% at 25 °C, 100 ± 0.1 % at 25 °C, certified to 1000 µS/cm)

Conductivity	Conductivity	Conductivity
0.720	0.713 (Reference)	1.000 (Reference)
0.719	0.713	0.719

The stated measurement uncertainty is given as the standard uncertainty multiplied with a coverage factor of 2 (i.e. 95% confidence level) and is stated in accordance with ISO 9000. The above uncertainty has been calculated at a confidence level of 95%.

The certificate is issued to the order.

The calibration has been performed under the following conditions: 25 °C, 50% RH.

Part of this calibration certificate can only be reproduced with the written consent of DFM.

Details of use of the certificate to be published in accordance with the manufacturer's instructions.

This certificate is consistent with the conditions that are included in Appendix 2 of the BIPM Guide to the Units, under the heading "4.1.2.4. Calibration certificates for the units of mass, length, time, and temperature".

Details of use of the certificate to be published in accordance with the manufacturer's instructions.

Details of use of the certificate to be published in accordance with the manufacturer's instructions.

2008-08-13

 Dr. Michael Maynard

DFM is accredited to ISO 9001:2008 for Calibration and Maintenance of Instruments
 DFM is accredited to ISO 9001:2008 for Calibration and Maintenance of Instruments

Certificate no. 02662
 Page 1 of 2
 Date: 2008-08-13

Recall:
 The solution was supplied by the client. Solution samples were provided in gas bottles containing 300 ml. The bottles were closed with a white cap and seal.

The conductivity has been determined from a measurement of the resistance of the solution in the conductivity bridge from 200 to 1000 Hz and from the calculated cell constant of the cell used. The cell was temperature controlled by submersion in an oil bath with stability better than ± 0.01 °C. The value for the conductivity has been referred to the given temperature (at 25 °C) by a correction based on the measurement of the temperature coefficient of the solution, determined at 25 °C.

The conductivity cell was last calibrated 2007-10-31 (DMF certificate 023024).
 The temperature bridge was last calibrated 2008-01-13 (DMF certificate 028002).
 The measurement was last calibrated 2008-01-28 (DMF certificate 03041).
 No conductivity contribution for CO₂ solubility of the solution has been taken into account.
 The average CO₂ partial pressure was measured near the measurement position and is given with the result.

DMF participates in the international collaboration under the Metro Convention in the CCEM (International Comparison of CO₂ solubility). This includes participation in international comparison of results and of uncertainty.

International comparison results and approved measurement capabilities are available at <http://bipm.org/>.

DMF (Data Measurement Facility)
 The Data Measurement Facility (DMF) is a not-for-profit research centre based at a campus of the University of Cambridge. The DMF provides a service to the scientific community in the UK and abroad in the field of metrology. The DMF is a member of the BIPM (International Bureau of Weights and Measures) and is a participant in the BIPM's international comparisons of results and of uncertainty. The DMF is also a member of the CCEM (International Comparison of CO₂ solubility). The DMF is a participant in the CCEM's international comparison of results and of uncertainty. The DMF is also a participant in the CCEM's international comparison of results and of uncertainty. The DMF is also a participant in the CCEM's international comparison of results and of uncertainty. The DMF is also a participant in the CCEM's international comparison of results and of uncertainty. The DMF is also a participant in the CCEM's international comparison of results and of uncertainty.

Value at 25°C	Accuracy	Stability (Months)	Certified by	Package	Order No.
1.3 µS/cm	± 1%	12	DFM	Glass bottle 300 ml	238 973
5 µS/cm	± 1%	36	DFM	Glass bottle 300 ml	238 926
15 µS/cm	± 1%	36	DFM	Glass bottle 300 ml	238 927
84 µS/cm	± 1%	18	DFM	Calpack bottle 500 ml	238 984
100 µS/cm	± 1%	36	DFM	Glass bottle 300 ml	238 934
147 µS/cm	± 1%	18	DFM	Calpack bottle 500 ml	238 985
1413 µS/cm	± 1%	36	DFM	Glass bottle 300 ml	238 928
1413 µS/cm	± 1%	18	DFM	Calpack bottle 500 ml	238 986
12880 µS/cm	± 1%	18	DFM	Calpack bottle 500 ml	238 988

Calibration Solutions

Potassium Chloride Solutions (Conductivity in mS/cm)

Temperature	Concentration ¹⁾		
°C	0.01 mol/l	0.1 mol/l	1 mol/l
0	0.776	7.15	65.41
5	0.896	8.22	74.14
10	1.020	9.33	83.19
15	1.147	10.48	92.52
16	1.173	10.72	94.41
17	1.199	10.95	96.31
18	1.225	11.19	98.22
19	1.251	11.43	100.14
20	1.278	11.67	102.07
21	1.305	11.91	104.00
22	1.332	12.15	105.94
23	1.359	12.39	107.89
24	1.386	12.64	109.84
25	1.413	12.88	111.80
26	1.441	13.13	113.77
27	1.468	13.37	115.74
28	1.496	13.62	
29	1.524	13.87	
30	1.552	14.12	
31	1.581	14.37	
32	1.609	14.62	
33	1.638	14.88	
34	1.667	15.13	
35	1.696	15.39	
36		15.64	

1) Data source: K. H. Hellwege (Editor), H. Landolt, R. Börnstein: Zahlenwerte und Funktionen ..., volume 2, part. volume 6


2) Data source: Test solutions calculated according to DIN IEC 746-3

Sodium Chloride Solutions (Conductivity in mS/cm)

Temperature	Concentration		
°C	0.01 mol/l ⁽²⁾	0.1 mol/l ⁽²⁾	Saturated ⁽¹⁾
0	0.631	5.786	134.5
1	0.651	5.965	138.6
2	0.671	6.145	142.7
3	0.692	6.327	146.9
4	0.712	6.510	151.2
5	0.733	6.695	155.5
6	0.754	6.881	159.9
7	0.775	7.068	164.3
8	0.796	7.257	168.8
9	0.818	7.447	173.4
10	0.839	7.638	177.9
11	0.861	7.831	182.6
12	0.883	8.025	187.2
13	0.905	8.221	191.9
14	0.927	8.418	196.7
15	0.950	8.617	201.5
16	0.972	8.816	206.3
17	0.995	9.018	211.2
18	1.018	9.221	216.1
19	1.041	9.425	221.0
20	1.064	9.631	226.0
21	1.087	9.838	231.0
22	1.111	10.047	236.1
23	1.135	10.258	241.1
24	1.159	10.469	246.2
25	1.183	10.683	251.3
26	1.207	10.898	256.5
27	1.232	11.114	261.6
28	1.256	11.332	266.9
29	1.281	11.552	272.1
30	1.306	11.773	277.4
31	1.331	11.995	282.7
32	1.357	12.220	288.0
33	1.382	12.445	293.3
34	1.408	12.673	298.7
35	1.434	12.902	304.1
36	1.460	13.132	309.5


Concentration Curves

-01- Sodium chloride solution NaCl



Conductivity versus substance concentration and process temperature for sodium chloride solution (NaCl)

-02- Hydrochloric acid HCl




Conductivity versus substance concentration and process temperature for hydrochloric acid (HCl)

Source: Haase/Sauermann/Dücker; Z. phys. Chem. New Edition, Vol. 47 (1965)


Concentration Curves

-03- Sodium hydroxide solution NaOH



Conductivity versus substance concentration and process temperature for sodium hydroxide solution (NaOH)

-04- Sulfuric acid H_2SO_4




Conductivity versus substance concentration and process temperature for sulfuric acid (H_2SO_4)

Source: Darling; Journal of Chemical and Engineering Data; Vol.9 No. 3, July 1964

Concentration Curves

-05- Nitric acid HNO_3



Conductivity versus substance concentration and process temperature for nitric acid (HNO_3)

Source: Haase/Sauermann/Dücker; Z. phys. Chem. New Edition, Vol. 47 (1965)

Conductance	Conductance $G [S] = 1 / R [\Omega]$
Conductivity	Conductivity $\chi [S/cm] = G [S] \cdot c [1/cm]$
Conductivity sensor	Either 2- or 4-electrode sensors can be connected. The cell constant of the sensor in use must be entered or be determined using a calibration solution taking account of the temperature.
Passcode	Preset four-digit number to select certain functions.
Sensocheck	Sensocheck monitors the sensor and its wiring. The resulting information is indicated by the Sensoface smileys. Sensocheck can be switched off.
Sensoface	Provides information on the sensor condition. Significant sensor polarization effects or an excessive cable capacitance are indicated..
Temperature coefficient	With temperature compensation activated, the measured value is calculated to the value at the reference temperature (25 °C) using the temperature coefficient.
Temperature compensation	Calculates the measured conductivity value for a reference temperature.

Index

22 mA signal for error message, output 1 41
22 mA signal for error message, output 2 47

A

Accessories 77
Alarm contact, configuration 51
Alarm contact, specifications 80
Alarm settings, configuration 50
Alarm settings, description 25
Assembly 10

C

Calibration 60
 by entry of cell constant 62
 Product calibration 66
 Temp probe adjustment 68
 with calibration solution 64
Calibration solutions 83
Cal Info 69
"Clean" contact 54
Concentration curves 86
 Hydrochloric acid HCl 87
 Nitric acid HNO₃ 90
 Sodium chloride solution NaCl 86
 Sodium hydroxide solution NaOH 88
 Sulfuric acid H₂SO₄ 89
Concentration measurement 35
Configuration 26
 Menu structure 27
Configuration: Alarm settings 50
Configuration: Limit function 52
Configuration: Output 1
 Concentration measurement 34
 Output current during Error and HOLD 40
 Output current range 36
 Output filter 38
 Process variable 32
 Sensor type 30

- Configuration: Output 2
 - Output current 42
 - Output current during HOLD 46
 - Output filter 44
 - Temperature 42
 - Temperature error 46
- Configuration: Rinsing probes 54
- Configuration: Temperature compensation 48
- Connection 14
- Current start / end 37, 43

D

- Default settings 56
- Device self-test 23
- Diagnostics functions 69
 - Display of calibration data 69
 - Display of last error message 69
 - Display of output currents 69
 - Sensor monitor 69
 - Specifying the output current 70
- Display 21
- Disposal 2
- Documentation 8

E

- Error Info 69
- Error messages 71
 - Display of last error message 69

F

- Factory settings of parameters 56
- Front panel 20

G

- Glossary 91

Index

H

- HAMILTON conductivity standard 83
- HOLD mode 24
 - Output signal during HOLD, output 1 41
 - Output signal during HOLD, output 2 47
- Hysteresis 53

I

- Installation 14
- Intended use 7

K

- Keypad 22

M

- Measurement 68
- Mounting plan 11

O

- Operating states 73
- Output filter, output 1 38
- Output filter, output 2 44
- Overview of configuration steps 28
- Overview of transmitter 9

P

- Package contents 10
- Panel mounting 13
- Panel-mount kit 13
- Parameters
 - Factory settings 56
 - Individual settings 58
- Passcodes 96
- Pipe mounting 12
- Pipe-mount kit 12
- Product calibration 66
- Product line 77
- Protective hood 12
- Protective wiring 18

Q

Quickstart guides 8

R

Relay, cleaning function 54
Relay contacts, protective wiring 18
Relay, limit function 52
Rinse contact 54
Rinsing interval 55
Rinsing probes 54
Rinsing system 55

S

Safety information 5
Self-test 23
Sensocheck 74
 Configuration 51
Sensoface 74
Sensor connection 17
Sensor monitor 69
Specifications 78

T

Technical data 78
Technical terms 91
Temperature compensation 48
Temperature measurement, configuration 43
Temperature probe adjustment 68
Terminal assignments 14
Time constant of output filter, output 1 39
Time constant of output filter, output 2 45

U

User interface 20

V

VP cable connection 17

W

Warranty 2
Wiring of Conducell sensors 17

Passcodes

Calibration

Key + passcode	Menu item	Page
cal + 0000	CAL info (display of cell constant)	69
cal + 0110	Calibration (with standard solution)	64
cal + 1100	Cell constant adjustment	62
cal + 1105	Product calibration	66
cal + 1015	Temp probe adjustment	68

Configuration

Key + passcode	Menu item	Page
conf + 0000	Error info (display of last error, erase)	69
conf + 1200	Configuration	26
conf + 2222	Sensor monitor (resistance, temp)	69
conf + 5555	Current source 1 (specify output current)	70
conf + 5556	Current source 2 (specify output current)	70
conf + ▶ + 4321	Factory setting	56