DE

Oxysens™ Bedienungsanleitung

Einleitung

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Oxysens™ Operating Instructions

Introduction

These operating instructions refer to Hamilton's Oxysens maintenance-free oxygen sensors.

Ref	Name	Distinguishing features
237150	Oxysens	fixed cable, NTC 22 kOhm

Hamilton sensors are quality products manufactured with the latest technology. You can only attain the highest accuracy and a maximum standing time by exactly following the instructions below.

These instructions should be read, understood and followed by all staff using the device. Hamilton can assume no responsibility for damage and operational disruptions arising from failure to observe these instructions.

Liability

The Liability of Hamilton Bonaduz AG is detailed in the document «General Terms and Conditions of Sale and Delivery (GTS)», chapter 12.

Hamilton is expressly not liable for direct or indirect losses arising from use of the sensors. It must in particular be insured in this conjunction that malfunctions can occur on account of the inherently limited useful life of sensors contingent upon their relevant applications. The user is responsible for the calibration, maintenance and regular replacement of the sensors. In the case of critical sensor applications, Hamilton recommends using back-up measuring points in order to avoid consequential damages. The user is responsible for taking suitable precautions in the event of a sensor failure.

Intended use

Oxysens are maintenance-free sensors for the measurement of dissolved oxygen in aqueous media. The sensor is either screwed into the process using the PG13,5 thread, or immersed in a tub using the immersing set. Typical applications include water resources management and fish farming.

The Oxysens sensor has a built-in temperature sensor (NTC 22 kOhm). This temperature sensor is to be used only to compensate the oxygen signal and not to control process temperature.

If the sensor is used in potentially explosive atmospheres, refer to the section «Instructions for use in potentially explosive atmospheres».

Safety Precautions

This sensor is to be used only as intended and is to be kept in a condition that ensures complete safety. The specifications given in the section «Technical data» as regards temperature, pressure etc. may under no circumstances be exceeded.

Inappropriate use or misuse can be dangerous.

Assembly and maintenance may only be done by trained personnel. Since the inner part of the sensor is glass, handle it with care. The top of the sensor can be damaged by impacts.

If mechanical damage to the sensor occurs, electrolyte may leak out. Contact of the electrolyte with skin, eyes or mucous membrane is to be avoided (see Safety Data Sheet Ref 608914)

Ensure that when screwed into the process the PG13,5 thread and the O-ring incur no damage. O-rings are consumable items and must be replaced regularly – at least once a year. Even where all necessary safety measures have been complied with, there is still a further risk involving leaks or mechanical damage to the armature. Where there are seals or screws, gases or liquids could leak out undetected.

Before dismounting the sensor, make sure that there is no over-pressure in the process.

Every product for shipment or sent back for repair must be decontaminated. If working with hazardous liquids observe and carry out the maintenance procedures, paying particular attention to cleaning and decontamination. If the product becomes contaminated with biohazardous, radioactive or chemical material, it should be cleaned.

Initial operation

During unpacking please check for possible mechanical defects. In case of a complaint please return the sensor in its original packing to Hamilton.

Preparation for the measurement

1. Connect to a suitable amplifier. Follow operating instructions for the amplifier. Functions of Oxysens cable cores:
 - colorless/center: cathode
 - brown: anode
 - yellow and blue: temperature sensor
 - outer shield: no connection to sample.
2. Remove watering cap (1). Keep sensor in air with membrane (2) pointing downwards. Pat membrane dry (see Figure 1).
3. Turn on amplifier and wait for about 15 min. until the measuring value is stabilized.
4. Run calibration on the amplifier, e.g. adjust to 100%.
5. Screw the sensor into the housing with the PG thread. Ensure that the torque of 1.5 Nm will not be exceeded to avoid damaging of the O-ring. A little grease on sealing facilitates assembly. Keep sensor membrane (2) pointing downwards. The sensor is now ready for measurement tasks.

Air/gas bubbles may stick to the sensitive area of the sensor. As a consequence, the measurement value might be wrong. In case of a shutdown of the power supply (230V) the measurement value could be wrong.

Oxysens sensors are maintenance-free. Replace the sensor if it can no longer be calibrated.

Cleaning

If working with hazardous liquids observe and carry out the maintenance procedures, paying particular attention to cleaning and decontamination.

Cleaning, assembly and maintenance should be performed by personnel trained in such work. Do not use any abrasive tissues or cleaning materials and do not use any cleaning chemicals other then described above. Before removing the sensor from the measuring setup, always make sure that the setup is pressureless and cold and that no process medium can be accidentally spilled. When removing and cleaning the sensor, it is recommended to wear safety glasses and protective gloves.

Accessories

Immersing Set: Hamilton Ref 237158, immerses the sensor and increases mechanical protection.

Picture B: Oxysens mounted in Immersing Set.

Instructions for use in potentially explosive atmospheres

ATEX / UKEX / IECEx marking:

Gas: CE 0035 / UKCA 2571
Ex II 1/2 G Ex ia IIC T4/T5/T6 Ga/Gb

Dust: CE 0035 / UKCA 2571
Ex II 1/2 D Ex ia IIIC T x °C Da/Db

Manufacturer:
Hamilton Bonaduz AG, 7402 Bonaduz, Switzerland

EU-Type Examination Certificate: TÜV 03 ATEX 7005 X
UK Type Examination Certificate: TÜV 21 UKEX 7047 X
IECEx Certificate of Conformity: IECEx TUR 14.0001 X

EU-Type Examination Certificate, UK Type Examination Certificate, and IECEx Certificate of Conformity can be downloaded from www.hamiltoncompany.com.

Conditions described in the ATEX EU-Type Examination Certificate, UK Type Examination Certificate, or IECEx Certificate of Conformity must be adhered to.

ATTENTION! In case a gas atmosphere and a dust atmosphere are or could be present at the same time, the risk of explosion must be examined carefully and special precautions may be necessary.

Coupling situation

a) The operator of equipment in potentially explosive atmospheres is responsible for ensuring that all components of the system are certified for that area classification and are compatible with each other.

b) Sensors, transmitter and other required equipment are to be set up within a single potential equalizing system.

c) O-rings having a sealing function between Zone 0 and 1, 20 and 21 respectively must be replaced at each dismantling of the sensor.

d) When the sensor is assembled, the attached ATEX / UKEX / IECEx sticker is to be attached to the cable in an easily visible place, as near as possible to the sensor itself. This sticker should not be removed at any time.

e) In gas atmospheres, the ex approval is not dependent on any conditions regarding mounting of the sensor. In dust atmospheres, however, special restrictions may apply when using armatures made of plastic material.

ATTENTION! Inappropriate fitting or dismantling of the armature or the sensor may lead to unintended release of an explosive atmosphere.

Connection to transmitter

Oxysens is suitable for connection to an intrinsically safe electrical circuit with protection level ia. The operator of the equipment must ensure that the allowable electrical values for the sensor all exceed those for the transmitter. None of the values cited for electrical power, voltage and current may be exceeded in total (measurement and temperature circuit together). The electrical values for Oxysens are as follows:

U = 24 V; I = 173 mA; P: Function of process temperature

For gas atmospheres:
The maximum allowable electric power P of the transmitter (measurement and temperature circuit together) depends on the desired process temperature.

T-class	P ≤ 30 mW	P ≤ 60 mW	P ≤ 90 mW
T4	60 °C	60 °C	60 °C
T5	60 °C	60 °C	60 °C
T6	60 °C	57 °C	46 °C

For dust atmospheres:
Table for the calculation of the maximum surface temperature «x» of the sensor as a function of the maximum electrical power of the transmitter «P_i» and ambient/process temperature «T_a». «x» must be smaller than the ignition temperatures of the dust involved.

Pi ≤ 30 mW	Pi ≤ 60 mW	Pi ≤ 90 mW
x=Ta + 12°C	x=Ta+23°C	x=Ta+34°C

The temperature limit (60°C) must always be respected.

Technical data

Electrode system	Silver platinum combination
Electrolyte	Oxylyte
Membrane / Cap	Delivered with Optiflow
Temperature Sensor	NTC 22 kOhm
Polarization Voltage and Time	-670 ± 50 mV, ≥ 2 hours
Current in air (at 25 °C)	40 - 80 nA
Response Time t98%	30 to 60 s at 25 °C, from air to nitrogen
Operating Temperature Range	0 to 60 °C
Storage temperature	-10 to 60 °C, with water containing watering cap
Pressure Range	0 to 4 bar inserted; max. 0.5 bar totally immersed
Shaft diameter	12 mm
Process Connection	PG13,5
Wetted Parts	Stainless steel 1.4435, Silicone, PEEK, FPM
Regeneration	not required
Temperature response	ca. 3.1%/K
Required Flow	≥ 0.03 m/s
Flow dependence	< 5% at 25 °C
Oxygen consumption	Ca. 20 ng/h in air at 25°C
Residual current	< 0.5% of current in air
Zero shift	< 0.5% of current in air every 2 months at 25°C in water under stable conditions
Sensitivity shift	< 5% every 2 months at 25°C in water under stable conditions


Figure 1

Disposal

The design of Hamilton sensors minimizes environmental impact. According to the EU directive 2012/19/EU the Hamilton sensors should be disposed as waste of electrical and electronic equipment, and not in municipal waste, or it can be sent back to Hamilton for disposal.

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© 2023 Hamilton Bonaduz AG. All rights reserved.

REF 610668/09 — 02/2023

HAMILTON

Web: www.hamiltoncompany.com
USA: 800-648-5950
Europe: +41-58-610-10-10

Hamilton Americas & Pacific Rim
4970 Energy Way
Reno, Nevada 89502 USA
Tel: +1-775-858-3000
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41-58-610-10-10
contact.pa.ch@hamilton.ch

To find a representative in your area, please visit www.hamiltoncompany.com.

This guide may be available in other languages.
Visit www.hamiltoncompany.com for more information.