




This drawing is the property of HAMILTON COMPANY and is furnished in confidence for a specific purpose. Disclosure or use of its contents for any other purpose is prohibited unless authorized by HAMILTON COMPANY.

Model Nr.:	242585	Model Rev.:	08-7-Wird bearbeitet
Object type:	ASSEM	Drawing Rev.:	04-4-Freigegeben



Pipe connection	L	ϕd
TC25 1/2"	19,3	9,4
TC25 3/8"	19,3	7,7
TC25 1/4"	19,3	4,6
Swagelok 1/2"	42,1	5,9
Swagelok 3/8"	39,3	5,9
Swagelok 1/4"	37,8	4,8
Swagelok 10	39,3	5,8
Swagelok 6	37,5	4,8

Pipe connection on both sides are identical



Information in this document is valid only for the revision indicated and may be subject to change without notice!

Dimensions without tolerances are for information purpose only!

Title: FlowCell	Material:	Mass:		
Created / Changed: Antl Christian 23929	Case-/Change Nr.: 11011173	Unit: [mm]	Format: A3	Rev.: 04
Additional Title:	Approved:	Release date:	Scale: 1/2	
Document number: SZC242585				