# Declaration of Quality

# EasyFerm Plus PHI VP 225 Pt100



#### **Product Identification**

Ref	238633-1141	
SN	37618	
Lot	111200384	

### **Factory Test Conditions**

Tests were performed according to	606281	*****
Buffer Temperature	24.6 °C	****

## **Factory Test Results**

	Test Specification	Result
Reading in pH7 (zero point)	0 mV ± 20 mV	6 mV
Reading in pH4		184 mV
Slope (at 25°C)	> 97 %	passed
Quality Test successfully passed	2021-02-08	



This declaration is generated automatically and valid without signature.



