

MATERIAL SPECIFICATION

Declaration of Compliance according to EN 10204-2.1

O-Rings at PG 13,5

REFs 10181562, 10191075, 10212283

Hamilton AG, Bonaduz, Switzerland

hereby declares that the wetted parts fulfill the following material specifications:

Wetted Parts	Material	Regulatory/Specification
O-ring(s) at PG 13,5	EPDM	FDA 21 CFR 177.2600 EC 1935/2004 ISO 10993-5 USP <87> USP <88> class VI (121 °C)
O-ring(s) at PG 13,5	FKM	FDA 21 CFR 177.2600 EC 1935/2004 ISO 10993-5 USP <87> USP <88> class VI (121 °C)
O-ring(s) at PG 13,5	VMQ (Silicone elastomer)	FDA 21 CFR 177.2600 EC 1935/2004 USP <87> USP <88> class VI (121 °C)

Hamilton AG

Arquint
Philipp 170



Digital unterschrieben
von Arquint Philipp 170
Datum: 2026.09.30
11:11:33 +02'00'

Dr. Philipp Arquint

Vice President Process Analytics

Pianegonda
Peter 2160



Digitally signed by
Pianegonda Peter 2160
Date: 2026.10.02 07:43:07
+02'00'

Peter Pianegonda

Senior Quality Manager
PA QA RA



To find a representative in your area, please visit:
hamiltoncompany.com/contact

Hamilton Americas & Pacific Rim
4970 Energy Way
Reno, Nevada 89502 USA
Tel: +1-775-858-3000
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41-58-610-10-10
contact.pa.ch@hamilton.ch