




This drawing is the property of HAMILTON COMPANY and is furnished in confidence for a specific purpose. Disclosure or use of its contents for any other purpose is prohibited unless authorized by HAMILTON COMPANY.

Model Nr.:	243900	Model Rev.:	00-3-Freigegeben
Object type:	ASSEM	Drawing Rev.:	00-2-Freigegeben

All required tube configurations possible



Pipe connection	L	ϕd
TC25 1/4"	16	4,57
TC25 3/8"	16	7,75
TC25 1/2"	16	9,4
TC25 3/4"	16	15,7
TC50 1"	16	22,1
TC50 1,5"	21	34,8
TC64 2"	24	47,5



Information in this document is valid only for the revision indicated and may be subject to change without notice!

Dimensions without tolerances are for information purpose only!

Title: Flow Cell UR		Material:		Mass:			
Additional Title:		Created / Changed: Antl Christian 23929		Case-/Change Nr.: 11012516		Unit: [mm]	Format: A3
		Approved: -		Release date: -		Scale: 1/2	Sheet: 1 / 1
		HAMILTON Bonaduz AG Via Crusch 8 CH-7402 Bonaduz		Document number: SZC243900		Rev.: 00	